

# Anti-Hairy Cell Leukemia antibody [DBA.44], prediluted ab857

## 1 图像

### 概述

产品名称	Anti-Hairy Cell Leukemia抗体[DBA.44], prediluted
描述	小鼠单克隆抗体[DBA.44] to Hairy Cell Leukemia, prediluted
宿主	Mouse
特异性	This antibody reacts with over 97% of hairy cell leukaemias as well as about 35% of high grade B cell lymphomas. It reacts with B lymphocytes, cells of mantle zone and immunoblasts outside the lymphoid tissues.
经测试应用	适用于: IHC-P
种属反应性	与反应: Human
免疫原	Tissue, cells or virus corresponding to Hairy Cell Leukemia. BALB/C mice were injected with the Deau cell line established from a large cell lymphoma.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C. Do Not Freeze.
存储溶液	pH: 7.3 Preservative: 0.05% Sodium azide
克隆	Inert stabilizer 单克隆
克隆编号	DBA.44
骨髓瘤	unknown
同种型	IgM

轻链类型

kappa

应用

The Abpromise guarantee

Abpromise™承诺保证使用ab857于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。


应用	Ab评论	说明
IHC-P		Use at an assay dependent concentration.

靶标

细胞定位

Nuclear

图片



Immunohistochemistry (Formalin fixed paraffin-embedded sections) analysis of human tonsil tissue labelling Hairy Cell Leukemia with ab857.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Hairy Cell Leukemia antibody [DBA.44], prediluted (ab857)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors