

Anti-Hantavirus Puumala strain antibody [1D11] ab116453

概述

产品名称	Anti-Hantavirus Puumala strain抗体[1D11]
描述	小鼠单克隆抗体[1D11] to Hantavirus Puumala strain
宿主	Mouse
经测试应用	适用于: WB, ELISA
种属反应性	与反应: Puumala hantavirus
免疫原	Recombinant Hantavirus Puumala strain.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	pH: 7.40 Constituent: 99% PBS
纯度	Protein A purified
纯化说明	Purity is > 95% by SDS-PAGE.
克隆	单克隆
克隆编号	1D11
同种型	IgG2b
轻链类型	kappa

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用ab116453于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 1 kDa.
ELISA		Use at an assay dependent dilution.

靶标

相关性 Hantaviruses are rodent-borne agents that cause severe human diseases. The viruses constitute a group of antigenetically and genetically related envelope viruses with a tripartite, negative-sense RNA genome. In northern, western and central Europe the Puumala virus (PUUV) is a predominant serotype. It is the etiological agent of nephropathia epidemica, a mild form of hemorrhagic fever with renal syndrome.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors