

HRP Anti-Apolipoprotein B antibody ab27622

3 References

概述

产品名称	HRP Anti-Apolipoprotein B抗体
描述	HRP山羊多克隆抗体to Apolipoprotein B
宿主	Goat
偶联物	HRP
经测试应用	适用于: ELISA, WB
种属反应性	与反应: Human
免疫原	Low Density Lipoproteins containing ApoB-100
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	<p>pH: 6.50</p> <p>Preservative: 0.01% Thimerosal (merthiolate)</p> <p>Constituents: 0.134% PBS, 0.58% Sodium chloride, 1% BSA</p>
纯度	Affinity purified
纯化说明	Purified by human apo B-100-Sepharose™ affinity column.
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用ab27622于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		Use at an assay dependent concentration.
WB		Use at an assay dependent concentration.

靶标

功能	Apolipoprotein B is a major protein constituent of chylomicrons (apo B-48), LDL (apo B-100) and VLDL (apo B-100). Apo B-100 functions as a recognition signal for the cellular binding and internalization of LDL particles by the apoB/E receptor.
疾病相关	Hypobetalipoproteinemia, familial, 1 Familial ligand-defective apolipoprotein B-100 Defects in APOB associated with defects in other genes (polygenic) can contribute to hypocholesterolemia.
序列相似性	Contains 1 vitellogenin domain.
翻译后修饰	Palmitoylated; structural requirement for proper assembly of the hydrophobic core of the lipoprotein particle.
细胞定位	Cytoplasm. Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors