

HRP Anti-GST antibody ab58626

★★★★★ [2 Abreviews](#) [18 References](#) [1 图像](#)

概述

产品名称	HRP Anti-GST抗体
描述	HRP山羊多克隆抗体to GST
宿主	Goat
偶联物	HRP
经测试应用	适用于: ELISA, WB, IHC-P
种属反应性	与反应: Species independent
免疫原	Full length protein corresponding to GST. Glutathione-S-Transferase [Schistosoma japonicum]
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 6.50 Preservative: 0.01% Gentamicin sulphate Constituents: 0.42% Tripotassium orthophosphate, 0.87% Sodium chloride, 1% BSA
纯度	Ion Exchange Chromatography
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab58626于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA	★★★★★ (1)	1/20000.
WB	★★★★★ (1)	1/2000 - 1/5000.
IHC-P		1/500 - 1/2000.

靶标

相关性 GST (Glutathione S-Transferase) is a 26kDa protein encoded by the parasitic helminth *Schistosoma japonicum* and widely used in the pGEX family of GST plasmid expression vectors as a fusion protein with foreign proteins.

图片



ab58626 staining GST in by ELISA. The sample was blocked in 1% milk for 1 hour at 18°C. The primary was diluted at 1/2000 in PBS-Tween20 buffer and incubated with sample for 1 hour at 18°C.

ELISA - HRP Anti-GST antibody (ab58626)

This image is courtesy of an anonymous Abreview

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors