

HRP Anti-IgE antibody [1A2] ab106493

1 图像

概述	
产品名称	HRP Anti-IgE 抗体[1A2]
描述	HRP 小鼠单克隆抗体[1A2] to IgE
宿主	Mouse
偶联物	HRP
特异性	ab106493 is specific for Human IgE. It reacts with IgE purified from blood but does not cross-react with Human anti-NP IgG1, IgG2, IgG3, IgG4, IgA and IgM
经测试应用	适用于: Sandwich ELISA
种属反应性	与反应: Human 不与反应: Mouse
免疫原	Mouse-Human chimeric anti-NP IgE
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能	
形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.01% Thimerosal (merthiolate) Constituent: PBS
纯度	Protein A purified
克隆	单克隆
克隆编号	1A2
同种型	IgG1

The Abpromise guarantee

Abpromise™承诺保证使用ab106493于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Sandwich ELISA		1/1000.

靶标


相关性

IgE is the class of antibodies produced in the lungs, skin, and mucous membranes. It may protect against parasite invasion, but it is a major factor in allergic reactions. The antigen-specific IgE interacts with mast cells and eosinophils, triggers the release of histamine, leukotrienes and other substances that lead to the itching, sneezing and congestion of allergies - and the life threatening respiratory distress of asthma and anaphylactic shock.

细胞定位

Cell Membrane

图片



ab106493 at 1/1000 dilution detecting IgE by sandwich ELISA. The capture antibody is an alternative mouse monoclonal (IgG1) raised against Mouse-Human chimeric anti-NP IgE.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors