

Anti-ICAM4 antibody ab112554

★★★★★ [1 Abreviews](#) [2 References](#) [1 图像](#)

概述

| | |
|-------|---|
| 产品名称 | Anti-ICAM4抗体 |
| 描述 | 兔多克隆抗体to ICAM4 |
| 宿主 | Rabbit |
| 经测试应用 | 适用于: IHC-P |
| 种属反应性 | 与反应: Human 预测可用于: Mouse, Rat  |
| 免疫原 | Recombinant fragment, corresponding to amino acids 43-210 of Human ICAM4 (BC000046). |
| 阳性对照 | Human fetal colon tissue. |
| 常规说明 | <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p> |

性能

| | |
|------|--|
| 形式 | Lyophilized:Add 200ul sterile distilled water |
| 存放说明 | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles. |
| 存储溶液 | pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 1% BSA, 98% PBS |
| 纯度 | Immunogen affinity purified |
| 克隆 | 多克隆 |
| 同种型 | IgG |

应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

| 应用 | Ab评论 | 说明 |
|-------|------|----------------|
| IHC-P | | 1/100 - 1/500. |

靶标

| | |
|-------|--|
| 功能 | ICAM proteins are ligands for the leukocyte adhesion protein LFA-1 (integrin alpha-L/beta-2). ICAM4 is also a ligand for alpha-4/beta-1 and alpha-V integrins. |
| 组织特异性 | Erythrocytes. |
| 序列相似性 | Belongs to the immunoglobulin superfamily. ICAM family. Contains 2 Ig-like C2-type (immunoglobulin-like) domains. |
| 翻译后修饰 | N- and O-glycosylated. |
| 细胞定位 | Secreted and Cell membrane. |

图片



ab112554, at a 1/100 dilution, staining ICAM4 in formalin-fixed, paraffin-embedded Human fetal colon tissue by Immunohistochemistry. Note: membrane staining.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-ICAM4 antibody (ab112554)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors