

### Anti-IGSF3 antibody ab97704

#### 1 图像

#### 概述

产品名称	Anti-IGSF3抗体
描述	兔多克隆抗体to IGSF3
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Human 预测可用于: Mouse, Rat 
免疫原	Recombinant fragment, corresponding to a sequence within amino acids 323-563 of Human IGSF3 (NP_001007238).
阳性对照	H1299 and Raji whole cell lysates; 293T, A431, HeLa, HepG2 and MOLT4 cell lysates.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	<p>pH: 7.00</p> <p>Preservative: 0.01% Thimerosal (merthiolate)</p> <p>Constituents: 1.21% Tris, 0.75% Glycine, 10% Glycerol (glycerin, glycerine)</p>
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

#### 应用

## The Abpromise guarantee

**Abpromise™** 承诺保证使用 ab97704 于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度；实际最佳的稀释度/浓度应由使用者检定。


应用	Ab评论	说明
WB		1/500 - 1/3000. Predicted molecular weight: 135 kDa.

## 靶标

**相关性** IGSF3, Immunoglobulin superfamily member 3, contains 8 Ig-like C2-type (immunoglobulin-like) domains. It is expressed in a wide range of tissues with high expression in placenta, kidney and lung.

**细胞定位** Membrane; Single-pass type I membrane protein

## 图片



**All lanes :** Anti-IGSF3 antibody (ab97704) at 1/1000 dilution

**Lane 1 :** H1299 whole cell lysate

**Lane 2 :** Raji whole cell lysate

Lysates/proteins at 30 µg per lane.

**Predicted band size:** 135 kDa

5% SDS-PAGE.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors