

# Anti-Interferon alpha 1 antibody [EPR18694-27] - BSA and Azide free (Capture) ab242952

**重组 RabMAb**

## 2 图像

### 概述

<b>产品名称</b>	Anti-Interferon alpha 1 抗体[EPR18694-27] - BSA and Azide free (Capture)
<b>描述</b>	兔单克隆抗体[EPR18694-27] to Interferon alpha 1 - BSA and Azide free (Capture)
<b>宿主</b>	Rabbit
<b>经测试应用</b>	<b>适用于:</b> Sandwich ELISA
<b>种属反应性</b>	<b>与反应:</b> Human
<b>免疫原</b>	Recombinant full length protein. This information is proprietary to Abcam and/or its suppliers.
<b>常规说明</b>	<p>ab242952 is a BSA and Azide Free antibody supplied in an unconjugated format and it is suitable for sandwich ELISAs to quantify Human Interferon alpha 1. The recommended pair for sandwich ELISA is:</p> <p style="padding-left: 20px;">Capture: ab242952, Human Interferon alpha 1 Capture Antibody (unconjugated)</p> <p style="padding-left: 20px;">Detector: <b>ab242984</b>, Human Interferon alpha 1 Detector Antibody (unconjugated)</p> <p>The reference range value is 2.34 - 150 pg/ml.</p> <p>Our <b>carrier-free</b> antibodies are typically supplied in a PBS-only formulation, purified and free of BSA, sodium azide and glycerol. The carrier-free buffer and high concentration allow for increased conjugation efficiency.</p> <p>This conjugation-ready format is designed for use with fluorochromes, metal isotopes, oligonucleotides, and enzymes, which makes them ideal for antibody labelling, functional and cell-based assays, flow-based assays (e.g. mass cytometry) and Multiplex Imaging applications.</p> <p>Use our <b>conjugation kits</b> for antibody conjugates that are ready-to-use in as little as 20 minutes with &lt;1 minute hands-on-time and 100% antibody recovery: available for fluorescent dyes, HRP, biotin and gold.</p> <p>The recommended antibody orientation is based on internal optimization for ELISA-based assays. Antibody orientation is assay dependent and needs to be optimized for each assay type. Please note that the range provided for this antibody is only an estimation based on the performance of the product using the recommended antibody pair. Performance of the antibody pair will depend on the specific characteristics of your assay. We guarantee the product works in sandwich ELISA, but we do not guarantee the sensitivity or dynamic range of the antibody in your assay.</p> <p>This product is a recombinant monoclonal antibody, which offers several advantages including:</p>

- High batch-to-batch consistency and reproducibility
- Improved sensitivity and specificity
- Long-term security of supply
- Animal-free production

For more information [see here](#).

Our RabMAb<sup>®</sup> technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to [RabMAb<sup>®</sup> patents](#).

## 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Constituent: 100% PBS
无载体	是
纯度	Protein A purified
克隆	单克隆
克隆编号	EPR18694-27
同种型	IgG

## 应用

**The Abpromise guarantee**      **Abpromise<sup>™</sup>承诺保证使用ab242952于以下的经测试应用**


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
<b>Sandwich ELISA</b>		Use at an assay dependent concentration. Can be paired for Sandwich ELISA with <b><u>Rabbit monoclonal [EPR18694-3] to Interferon alpha 1 - BSA and Azide free (Detector) (ab242984)</u></b> .

## 靶标

功能	Produced by macrophages, IFN-alpha have antiviral activities. Interferon stimulates the production of two enzymes: a protein kinase and an oligoadenylate synthetase.
序列相似性	Belongs to the alpha/beta interferon family.
细胞定位	Secreted.





## 图片



Representative standard curve from corresponding SimpleStep ELISA® Kit ([ab213479](#)).

Sandwich ELISA - Anti-Interferon alpha 1 antibody [EPR18694-27] - BSA and Azide free (Capture) (ab242952)

Why choose a recombinant antibody?

 <p><b>Research with confidence</b> Consistent and reproducible results</p>	 <p><b>Long-term and scalable supply</b> Recombinant technology</p>
 <p><b>Success from the first experiment</b> Confirmed specificity</p>	 <p><b>Ethical standards compliant</b> Animal-free production</p>

Anti-Interferon alpha 1 antibody [EPR18694-27] - BSA and Azide free (Capture) (ab242952)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors