

Anti-LYRM1 antibody ab106665

1 图像

概述

产品名称	Anti-LYRM1抗体
描述	兔多克隆抗体to LYRM1
宿主	Rabbit
特异性	Two isoforms of LYRM1 are known to exist; ab106665 will recognize both isoforms. ab106665 is predicted not to cross-react with other LYRM family members.
经测试应用	适用于: WB
种属反应性	与反应: Rat, Human
免疫原	A 17 amino acid synthetic peptide near the center of Human LYRM1 (EAW66844).
阳性对照	Human liver tissue lysate;
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at 4°C (stable for up to 12 months).
存储溶液	<p>pH: 7.2</p> <p>Preservative: 0.02% Sodium azide</p> <p>Constituent: PBS</p>
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用 ab106665 于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度；实际最佳的稀释度/浓度应由使用者检定。

应用	Ab 评论	说明
WB		Use a concentration of 1 - 2 µg/ml. Predicted molecular weight: 28 kDa.

靶标

功能	May promote cell proliferation and inhibition of apoptosis of preadipocytes.
组织特异性	High levels in adipose tissue.
序列相似性	Belongs to the complex I LYR family.
细胞定位	Nucleus.

图片



Lane 1 : Anti-LYRM1 antibody (ab106665) at 1 µg/ml

Lane 2 : Anti-LYRM1 antibody (ab106665) at 2 µg/ml

All lanes : Human liver tissue lysate

Lysates/proteins at 15 µg per lane.

Predicted band size: 28 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors