

Anti-MARCO antibody [EPR22944-64] ab239369

重组 RabMAb

8 图像

概述

产品名称	Anti-MARCO抗体[EPR22944-64]
描述	兔单克隆抗体[EPR22944-64] to MARCO
宿主	Rabbit
经测试应用	适用于: IP, IHC-Fr, IHC-P, ICC/IF, WB 不适用于: Flow Cyt
种属反应性	与反应: Mouse
免疫原	Synthetic peptide. This information is proprietary to Abcam and/or its suppliers.
阳性对照	WB: J774A.1 treated with Lipopolysaccharides whole cell lysate. Mouse spleen lysate. Mouse spleen lysate treated with Peptide: N-glycosidase F (PNGase F). IHC-P: Mouse spleen and liver tissue. IHC-Fr: Mouse spleen tissue. ICC/IF: J774A.1 cells treated with Lipopolysaccharides. IP: Mouse spleen lysate.
常规说明	This product is a recombinant monoclonal antibody, which offers several advantages including: <ul style="list-style-type: none"> - High batch-to-batch consistency and reproducibility - Improved sensitivity and specificity - Long-term security of supply - Animal-free production For more information see here . Our RabMAb [®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to RabMAb[®] patents .

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
存储溶液	pH: 7.2 Preservative: 0.01% Sodium azide Constituents: PBS, 40% Glycerol (glycerin, glycerine), 0.05% BSA
纯度	Protein A purified
克隆	单克隆

克隆编号

EPR22944-64

同种型

IgG

应用

The Abpromise guarantee

Abpromise™承诺保证使用ab239369于以下的经测试应用

“应用说明”部分下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IP		1/30.
IHC-Fr		1/500.
IHC-P		1/500. Perform heat mediated antigen retrieval with Tris/EDTA buffer pH 9.0 before commencing with IHC staining protocol.
ICC/IF		1/100.
WB		1/1000.

应用说明

Is unsuitable for Flow Cyt.

靶标

功能

Pattern recognition receptor (PRR). Binds Gram-positive and Gram-negative bacteria.

序列相似性


Contains 1 collagen-like domain.

Contains 1 SRCR domain.

细胞定位

Membrane.

图片




Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-MARCO antibody [EPR22944-64] (ab239369)

Immunohistochemical analysis of paraffin-embedded mouse liver tissue labeling MARCO with ab239369 at 1/500 dilution (1.11 µg/ml) followed by a ready to use Rabbit specific IHC polymer detection kit HRP/DAB (**ab209101**). Positive staining on Kupffer cells of mouse liver is observed. The section was incubated with ab239369 for 30 mins at room temperature. The immunostaining was performed on a Leica Biosystems BOND[®]RX instrument. Counterstained with hematoxylin.

Secondary antibody only control: Secondary antibody is a ready to use Rabbit specific IHC polymer detection kit HRP/DAB (**ab209101**).

Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0, epitope retrieval solution 2) for 20mins.



Western blot - Anti-MARCO antibody [EPR22944-64] (ab239369)

All lanes : Anti-MARCO antibody [EPR22944-64] (ab239369) at 1/1000 dilution

Lane 1 : Untreated J774A.1 (mouse reticulum cell sarcoma monocyte macrophage), whole cell lysate

Lane 2 : J774A.1 treated with 1ug/ml Lipopolysaccharides (LPS) for 24 hours, whole cell lysate

Lysates/proteins at 20 µg per lane.

Secondary


All lanes : Goat Anti-Rabbit IgG H&L (HRP) (**ab97051**) at 1/100000 dilution

Observed band size: 55-75 kDa

Blocking and dilution buffer: 5% NFDM/TBST .

Exposure time: 3 seconds.

The expression of MARCO/Macrophage receptor with collagenous structure is induced by Lipopolysaccharides (LPS) treatment (PMID: 9916718).



Western blot - Anti-MARCO antibody [EPR22944-64] (ab239369)

All lanes : Anti-MARCO antibody [EPR22944-64] (ab239369) at 1/1000 dilution

Lane 1 : Mouse spleen lysate

Lane 2 : Mouse spleen lysate treated with Peptide:N-glycosidase F (PNGase F)

Lysates/proteins at 20 µg per lane.

Secondary

All lanes : Goat Anti-Rabbit IgG H&L (HRP) (**ab97051**) at 1/100000 dilution


Observed band size: 53,75 kDa

Exposure time: 15 seconds

75 kDa (glycosylated form), 53 kDa (deglycosylated form).

Blocking and Dilution Buffer and concentration: 5% NFDm/TBST.

N-glycosylation is crucial for MARCO/Macrophage receptor with collagenous structure intercellular trafficking and surface targeting (PMID: 26892079).



Immunoprecipitation - Anti-MARCO antibody [EPR22944-64] (ab239369)

MARCO was immunoprecipitated from 0.35 mg mouse spleen lysate 10µg with ab239369 at 1/30 dilution (2µg in 0.35mg lysates).

Western blot was performed on the immunoprecipitate using ab239369 at 1/1000 dilution. VeriBlot for IP Detection Reagent (HRP) (**ab131366**) was used at 1/5000 dilution.


Lane 1: Mouse spleen lysate 10µg.

Lane 2: ab239369 IP in mouse spleen lysate.

Lane 3: Rabbit monoclonal IgG (**ab172730**) instead of ab239369 in mouse spleen lysate.

Blocking and dilution buffer and concentration: 5% NFDm/TBST.


Exposure time: 8 seconds.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-MARCO antibody [EPR22944-64] (ab239369)

Immunohistochemical analysis of paraffin-embedded mouse spleen tissue labeling MARCO with ab239369 at 1/500 dilution (1.11 µg/ml) followed by a ready to use Rabbit specific IHC polymer detection kit HRP/DAB (**ab209101**). Positive staining on macrophages of mouse spleen (PMID: 12874264, 15494482) is observed. The section was incubated with ab239369 for 30 mins at room temperature. The immunostaining was performed on a Leica Biosystems BOND[®] instrument. Counterstained with hematoxylin. Secondary antibody only control: Secondary antibody is a ready to use Rabbit specific IHC polymer detection kit HRP/DAB (**ab209101**).


Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0, epitope retrieval solution 2) for 20mins.



Immunohistochemistry (Frozen sections) - Anti-MARCO antibody [EPR22944-64] (ab239369)

4% PFA-fixed, 0.2% Triton X-100 permeabilized frozen section of mouse spleen tissue labeling MARCO using ab239369 at 1/500 dilution, followed by **ab150077** AlexaFluor[®]488 Goat anti-Rabbit secondary antibody at 1/1000 dilution (green). The nuclear counterstain is DAPI (blue). Positive staining on mouse splenic marginal zone (PMID: 24670793) is observed.

Heat mediated antigen retrieval using sodium citrate buffer (10mM citrate pH 6.0 + 0.05% Tween-20).



Immunocytochemistry/ Immunofluorescence - Anti-MARCO antibody [EPR22944-64] (ab239369)

Immunofluorescent analysis of 4% paraformaldehyde-fixed, 0.1% Triton X-100 permeabilized J744A.1 (Mouse reticulum cell sarcoma monocyte macrophage) cells labeling MARCO with ab239369 at 1/100 dilution, followed by **ab150077** AlexaFluor[®]488 Goat anti-Rabbit secondary at 1/1000 dilution (green). Confocal image showing cytoplasmic staining in J744A.1 cells treated with LPS (1 µg/ml) for 24h. The nuclear counterstain is DAPI (blue). Tubulin is detected with **ab195889** Anti-alpha Tubulin antibody [DM1A] - Microtubule Marker (Alexa Fluor[®] 594) at 1/200 dilution (red).

Secondary antibody only control: Used PBS instead of primary antibody, followed by **ab150077** AlexaFluor[®]488 Goat anti-Rabbit secondary at 1/1000 dilution.

Why choose a recombinant antibody?



Research with confidence
Consistent and reproducible results



Long-term and scalable supply
Recombinant technology



Success from the first experiment
Confirmed specificity



Ethical standards compliant
Animal-free production

Anti-MARCO antibody [EPR22944-64] (ab239369)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors