

# Anti-Measles nucleoprotein antibody ab23974

## 9 References

### 概述

---

产品名称	Anti-Measles nucleoprotein抗体
描述	兔多克隆抗体to Measles nucleoprotein
宿主	Rabbit
经测试应用	适用于: ELISA, WB
种属反应性	与反应: Measles virus
免疫原	Recombinant full length protein.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 性能

---

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Constituent: Whole serum
纯度	Whole antiserum
克隆	多克隆
同种型	IgG

### 应用

---

**The Abpromise guarantee**      **Abpromise™**承诺保证使用ab23974于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		Use at an assay dependent concentration.
WB		Use at an assay dependent concentration.

## 靶标

**相关性** Involved in the encapsidating of the genome in a ratio of 1 N per 6 ribonucleotides, resulting in the protection from nucleases. The nucleocapsid (NC) has an helical structure with either 12.35 or 11.64 N per turn, approximately 20nm in diameter, with a hollow central cavity approximately 5nm in diameter. The encapsidated genomic RNA is termed the NC and serves as template for transcription and replication. During replication, encapsidation by N is coupled to RNA synthesis and all replicative products are resistant to nucleases. N is released in the blood following lysis of measles infected cells, it interacts then with human FCGR2B on immune cells, inducing apoptosis and blocking inflammatory immune response. Ntail binds to a protein on human thymic epithelial cells, termed Nucleoprotein Receptor (NR), inducing growth arrest

**细胞定位** Virion. Host cytoplasm.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors