

Mouse IgG - Isotype Control ab37355

[65 References](#) [1 图像](#)

概述

产品名称	小鼠IgG -同型对照
特异性	This antibody has been selected as an isotype control as it has no known specificity.
经测试应用	适用于: ELISA, Flow Cyt, IHC
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
存储溶液	pH: 8.20 Constituent: 100% Borate buffered saline
纯度	Affinity purified
纯化说明	Purified from normal mouse serum.
克隆	多克隆
同种型	IgG

应用


The Abpromise guarantee **Abpromise™** 承诺保证使用ab37355于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		Use at an assay dependent concentration.

应用	Ab评论	说明
Flow Cyt		Use at an assay dependent concentration. ab37355 - Mouse polyclonal IgG, is suitable for use as an isotype control with this antibody.
IHC		Use at an assay dependent concentration.

图片



ELISA plate was coated with serially diluted Mouse IgG-UNLB (ab37355). Immunoglobulin was detected with Goat Anti-Mouse IgG, Human ads-BIOT, Goat Anti-Mouse IgA-BIOT, and Goat Anti-Mouse IgM, Human ads-BIOT followed by Streptavidin-HRP and quantified.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors