

Product datasheet

Mouse IgG2c, kappa monoclonal [18C8BC7AD10] - Isotype Control - BSA and Azide free ab221849

2 References

概述

产品名称	小鼠IgG2c, kappa单克隆抗体[18C8BC7AD10]-同型对照- BSA and Azide free
特异性	This Mouse IgG2c Isotype Control Antibody was raised in mouse against a yeast-specific protein making it unsuitable for yeast experiments.
经测试应用	适用于: Flow Cyt, ICC
常规说明	<p>This antibody clone is manufactured by Abcam.</p> <p>Ab221849 is a PBS only buffer formulated version of ab170191.</p> <p>Isotype controls are used to confirm that the primary antibody binding is specific and not a result of non-specific Fc receptor binding or other protein interactions. The isotype control antibody should match the primary antibody's host species, isotype, and possible conjugation.</p> <p>If you require this antibody in a particular buffer formulation or a particular conjugate for your experiments, please contact orders@abcam.com or you can find further information here.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 7.5 Constituent: 100% PBS
无载体	是
纯度	Ammonium Sulphate Precipitation

克隆	单克隆
克隆编号	18C8BC7AD10
同种型	IgG2c
轻链类型	kappa

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab221849于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		Use a concentration of 1 µg/ml.
ICC		Use a concentration of 1 µg/ml.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors