

Anti-NK-3R antibody ab124025

[2 Abreviews](#) [1 References](#) [2 图像](#)

概述

产品名称	Anti-NK-3R抗体
描述	兔多克隆抗体to NK-3R
宿主	Rabbit
特异性	ab124025 does not recognise NK1R or NK2R.
经测试应用	适用于: WB, IHC-P
种属反应性	与反应: Mouse, Human
免疫原	Synthetic peptide corresponding to Human NK-3R (internal sequence).
阳性对照	Human brain tissue; RAW 264.7 whole cell lysate
常规说明	<p>This product should be stored undiluted.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C long term.
存储溶液	Preservative: 0.02% Sodium azide Constituent: 99% PBS
纯度	Protein A purified
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee


Abpromise™承诺保证使用ab124025于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		Use a concentration of 0.5 - 1 µg/ml. Detects a band of approximately 46 kDa (predicted molecular weight: 52 kDa).
IHC-P		Use a concentration of 5 µg/ml. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

靶标	
功能	This is a receptor for the tachykinin neuropeptide neuromedin-K (neurokinin B). It is associated with G proteins that activate a phosphatidylinositol-calcium second messenger system. The rank order of affinity of this receptor to tachykinins is: neuromedin-K > substance K > substance P.
疾病相关	Hypogonadotropic hypogonadism 11 with or without anosmia
序列相似性	Belongs to the G-protein coupled receptor 1 family.
翻译后修饰	The anchoring of this receptor to the plasma membrane is probably mediated by the palmitoylation of a cysteine residue.
细胞定位	Cell membrane.

图片



Lane 1 : Anti-NK-3R antibody (ab124025) at 0.5 µg/ml

Lane 2 : Anti-NK-3R antibody (ab124025) at 2 µg/ml


All lanes : RAW 264.7 whole cell lysate

Predicted band size: 52 kDa

Observed band size: 46 kDa

Western blot - Anti-NK-3R antibody (ab124025)

2



ab124025, at 5µg/ml, staining NK-3R in paraffin-embedded Human brain tissue by Immunohistochemistry.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-NK-3R antibody (ab124025)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors