

Anti-PARK7/DJ1 antibody [EP2816Y] ab76241

敲除验证 RabMAb

★★★★★ 1 Abreviews 8 References 5 图像

概述

产品名称	Anti-PARK7/DJ1抗体[EP2816Y]
描述	兔单克隆抗体[EP2816Y] to PARK7/DJ1
宿主	Rabbit
经测试应用	适用于: WB, IHC-P, Flow Cyt (Intra) 不适用于: IP
种属反应性	与反应: Human 预测可用于: Mouse, Rat 
免疫原	Synthetic peptide. This information is proprietary to Abcam and/or its suppliers.
阳性对照	WB: TF-1, Jurkat, HEK293T, HAP1 and HeLa cell lysates; Human brain nuclear extract tissue lysate; Human brain tissue lysate. Flow Cyt (intra): Jurkat cells. IHC-P: Human brain tissue.
常规说明	<p>Our RabMAb<sup>®</sup> technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to <a href="#">RabMAb<sup>®</sup> patents</a>.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	pH: 7.20 Preservative: 0.05% Sodium azide Constituents: 0.1% BSA, 40% Glycerol (glycerin, glycerine), 9.85% Tris glycine, 50% Tissue culture supernatant
纯度	Tissue culture supernatant

克隆

单克隆

克隆编号

EP2816Y

同种型

IgG

应用

The Abpromise guarantee      **Abpromise™**承诺保证使用ab76241于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB	★★★★★ (1)	1/10000 - 1/20000. Detects a band of approximately 24 kDa (predicted molecular weight: 20 kDa).
IHC-P		1/100 - 1/250. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.
Flow Cyt (Intra)		1/50. <b>ab172730</b> - Rabbit monoclonal IgG, is suitable for use as an isotype control with this antibody.

应用说明      Is unsuitable for IP.

靶标

功能

Protects cells against oxidative stress and cell death. Plays a role in regulating expression or stability of the mitochondrial uncoupling proteins SLC25A14 and SLC25A27 in dopaminergic neurons of the substantia nigra pars compacta and attenuates the oxidative stress induced by calcium entry into the neurons via L-type channels during pacemaking. Eliminates hydrogen peroxide and protects cells against hydrogen peroxide-induced cell death. May act as an atypical peroxiredoxin-like peroxidase that scavenges hydrogen peroxide. Following removal of a C-terminal peptide, displays protease activity and enhanced cytoprotective action against oxidative stress-induced apoptosis. Stabilizes NFE2L2 by preventing its association with KEAP1 and its subsequent ubiquitination. Binds to OTUD7B and inhibits its deubiquitinating activity. Enhances RELA nuclear translocation. Binds to a number of mRNAs containing multiple copies of GG or CC motifs and partially inhibits their translation but dissociates following oxidative stress. Required for correct mitochondrial morphology and function and for autophagy of dysfunctional mitochondria. Regulates astrocyte inflammatory responses. Acts as a positive regulator of androgen receptor-dependent transcription. Prevents aggregation of SNCA. Plays a role in fertilization. Has no proteolytic activity. Has cell-growth promoting activity and transforming activity. May function as a redox-sensitive chaperone.

组织特异性


Highly expressed in pancreas, kidney, skeletal muscle, liver, testis and heart. Detected at slightly lower levels in placenta and brain. Detected in astrocytes, Sertoli cells, spermatogonia, spermatids and spermatozoa.

疾病相关

Defects in PARK7 are the cause of Parkinson disease type 7 (PARK7) [MIM:606324]. A neurodegenerative disorder characterized by resting tremor, postural tremor, bradykinesia, muscular rigidity, anxiety and psychotic episodes. PARK7 has onset before 40 years, slow progression and initial good response to levodopa. Some patients may show traits reminiscent of amyotrophic lateral sclerosis-parkinsonism/dementia complex (Guam disease).

序列相似性	Belongs to the peptidase C56 family.
翻译后修饰	Sumoylated on Lys-130 by PIAS2 or PIAS4; which is enhanced after ultraviolet irradiation and essential for cell-growth promoting activity and transforming activity. Cys-106 is easily oxidized to sulfinic acid. Undergoes cleavage of a C-terminal peptide and subsequent activation of protease activity in response to oxidative stress.
细胞定位	Cytoplasm. Nucleus. Mitochondrion. Under normal conditions, located predominantly in the cytoplasm and, to a lesser extent, in the nucleus and mitochondrion. Translocates to the mitochondrion and subsequently to the nucleus in response to oxidative stress and exerts an increased cytoprotective effect against oxidative damage. Detected in tau inclusions in brains from neurodegenerative disease patients.

图片



**All lanes :** Anti-PARK7/DJ1 antibody [EP2816Y] (ab76241) at 1/1000 dilution

- Lane 1 :** Wild-type HEK-293T (Human epithelial cell line from embryonic kidney transformed with large T antigen) whole cell lysate
- Lane 2 :** PARK7 knockout HEK-293T (Human epithelial cell line from embryonic kidney transformed with large T antigen) whole cell lysate
- Lane 3 :** HeLa (Human epithelial cell line from cervix adenocarcinoma) whole cell lysate
- Lane 4 :** Human brain nuclear fraction tissue lysate

Lysates/proteins at 20 µg per lane.

Secondary


**All lanes :** Goat anti-Rabbit IgG H&L (IRDye® 800CW) preadsorbed ([ab216773](#)) at 1/10000 dilution

**Predicted band size:** 20 kDa  
**Observed band size:** 24 kDa

**Lanes 1-4:** Merged signal (red and green). Green - ab76241 observed at 24 kDa. Red - loading control [ab8245](#) observed at 36 kDa.

ab76241 Anti-PARK7/DJ1 antibody [EP2816Y] was shown to specifically react with PARK7/DJ1 in wild-type HEK-293T cells. Loss of signal was observed when knockout cell line [ab266338](#) (knockout cell lysate [ab257016](#)) was used. Wild-type and

PARK7/DJ1 knockout samples were subjected to SDS-PAGE. ab76241 and Anti-GAPDH antibody [6C5] - Loading Control ([ab8245](#)) were incubated overnight at 4°C at 1 in 1000 dilution and 1 in 20000 dilution respectively. Blots were developed with Goat anti-Rabbit IgG H&L (IRDye® 800CW) preadsorbed ([ab216773](#)) and Goat anti-Mouse IgG H&L (IRDye® 680RD) preadsorbed ([ab216776](#)) secondary antibodies at 1 in 20000 dilution for 1 hour at room temperature before imaging.



**All lanes :** Anti-PARK7/DJ1 antibody [EP2816Y] (ab76241) at 1/10000 dilution

**Lane 1 :** Wild-type HAP1 cell lysate

**Lane 2 :** PARK/DJ1 knockout HAP1 cell lysate

**Lane 3 :** HeLa cell lysate

**Lane 4 :** Human brain tissue lysate

Lysates/proteins at 20 µg per lane.

### Secondary

**All lanes :** Goat anti-Rabbit IgG H&L (IRDye® 800CW) preadsorbed ([ab216773](#)) at 1/10000 dilution

**Predicted band size:** 20 kDa

**Observed band size:** 24 kDa

**Lane 1:** Wild-type HAP1 cell lysate (20 µg)

**Lane 2:** PARK/DJ1 knockout HAP1 cell lysate (20 µg)


**Lane 3:** HeLa cell lysate (20 µg)

**Lane 4:** Human brain tissue lysate (20 µg)

**Lanes 1 - 4:** Merged signal (red and green). Green - ab76241 observed at 24 kDa. Red - loading control, [ab8245](#), observed at 37 kDa.

ab76241 was shown to specifically react with PARK7/DJ1 in wild-type HAP1 cells. No band was observed when PARK7/DJ1 knockout samples were used. Wild-type and PARK7/DJ1 knockout samples were subjected to SDS-PAGE. ab76241 and [ab8245](#) (loading control to PARK7/DJ1) were both diluted 1/10,000 and incubated overnight at 4°C. Blots were developed with Goat anti-Rabbit IgG H&L (IRDye® 800CW) preadsorbed ([ab216773](#)) and Goat anti-Mouse IgG H&L (IRDye® 680RD) preadsorbed

([ab216776](#)) secondary antibodies at 1/10,000 dilution for 1hr at room temperature before imaging.



Western blot - Anti-PARK7/DJ1 antibody [EP2816Y] (ab76241)

**All lanes :** Anti-PARK7/DJ1 antibody [EP2816Y] (ab76241) at 1/20000 dilution

**Lane 1 :** TF-1 cell lysate

**Lane 2 :** Jurkat cell lysate

**Lane 3 :** HeLa cell lysate


Lysates/proteins at 10 µg per lane.

### Secondary

**All lanes :** HRP labelled goat anti-rabbit at 1/1000 dilution


**Predicted band size:** 20 kDa

**Observed band size:** 24 kDa



Flow Cytometry (Intracellular) - Anti-PARK7/DJ1 antibody [EP2816Y] (ab76241)

Overlay histogram showing Jurkat cells stained with ab76241 (red line). The cells were fixed with 80% methanol (5 min) and then permeabilized with 0.1% PBS-Tween for 20 min. The cells were then incubated in 1x PBS / 10% normal goat serum / 0.3M glycine to block non-specific protein-protein interactions followed by the antibody (ab76241, 1/100 dilution) for 30 min at 22°C. The secondary antibody used was Alexa Fluor® 488 goat anti-rabbit IgG (H+L) ([ab150077](#)) at 1/2000 dilution for 30 min at 22°C. Isotype control antibody (black line) was rabbit IgG (monoclonal) (1 µg/1x10<sup>6</sup> cells) used under the same conditions. Unlabelled sample (blue line) was also used as a control. Acquisition of >5,000 events were collected using a 20mW Argon ion laser (488nm) and 525/30 bandpass filter. This antibody gave a positive signal in Jurkat cells fixed with 4% paraformaldehyde (10 min)/permeabilized with 0.1% PBS-Tween for 20 min used under the same conditions.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-PARK7/DJ1 antibody [EP2816Y] (ab76241)

Immunohistochemical analysis of paraffin-embedded human brain tissue using ab76241 at 1/100 dilution.

Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors