




### Anti-PDAP1 antibody ab204340

#### 3 图像

#### 概述

产品名称	Anti-PDAP1抗体
描述	兔多克隆抗体to PDAP1
宿主	Rabbit
经测试应用	适用于: IHC-P, WB, ICC/IF
种属反应性	与反应: Human 预测可用于: Mouse, Rat 
免疫原	Recombinant fragment corresponding to Human PDAP1 aa 50-150. Database link: <a href="#">Q13442</a>  <a href="#">Run BLAST with</a>  <a href="#">Run BLAST with</a>
阳性对照	Human breast tissue; RT-4 and U-251 MG cell lysates; U-251 MG cells.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
存储溶液	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 59% PBS, 40% Glycerol (glycerin, glycerine)
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

## 应用

### The Abpromise guarantee

**Abpromise™** 承诺保证使用 ab204340 于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-P		1/200 - 1/500. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.
WB		1/250 - 1/500. Predicted molecular weight: 21 kDa.
ICC/IF		Use a concentration of 0.25 - 2 µg/ml.

## 靶标


### 功能

Enhances PDGFA-stimulated cell growth in fibroblasts, but inhibits the mitogenic effect of PDGFB.

### 序列相似性

Belongs to the PDAP1 family.

## 图片



Western blot - Anti-PDAP1 antibody (ab204340)


**All lanes :** Anti-PDAP1 antibody (ab204340) at 1/250 dilution

**Lane 1 :** RT-4 cell lysate

**Lane 2 :** U-251 MG cell lysate


Developed using the ECL technique.

**Predicted band size:** 21 kDa



Immunofluorescent analysis of PFA-fixed, Triton X-100 permeabilized U-251 MG cells labeling PDAP1 with ab204340 at 4 µg/ml (green).

Immunocytochemistry/ Immunofluorescence - Anti-PDAP1 antibody (ab204340)



Immunohistochemical analysis of paraffin-embedded Human breast tissue labeling PDAP1 with ab204340 at 1/200 dilution.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-PDAP1 antibody (ab204340)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

**Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors