

Anti-PDGF B antibody [MM0014-5F66] ab51869

★★★★★ [1 Abreviews](#) [6 References](#) [3 图像](#)

概述

产品名称	Anti-PDGF B抗体[MM0014-5F66]
描述	小鼠单克隆抗体[MM0014-5F66] to PDGF B
宿主	Mouse
经测试应用	适用于: WB, Flow Cyt, IHC-P
种属反应性	与反应: Human
免疫原	Recombinant full length protein (Human)
阳性对照	WB: Human placenta tissue lysate. IHC-P: Human placenta tissue. Flow Cyt: SHSY-5Y cells.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Constituent: PBS
纯度	Protein G purified
克隆	单克隆
克隆编号	MM0014-5F66
同种型	IgG1

应用

The Abpromise guarantee


Abpromise™ 承诺保证使用ab51869于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		1/1000. Predicted molecular weight: 27 kDa.
Flow Cyt		Use 1µg for 10 ⁶ cells. ab170190 - Mouse monoclonal IgG1, is suitable for use as an isotype control with this antibody.
IHC-P		1/50 - 1/200. Perform enzymatic antigen retrieval before commencing with IHC staining protocol.

靶标

图片




Western blot - Anti-PDGF B antibody [MM0014-5F66] (ab51869)

Anti-PDGF B antibody [MM0014-5F66] (ab51869) at 1/1000 dilution + Human placenta tissue lysate

Predicted band size: 27 kDa


Observed band size: 27 kDa



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-PDGF BB antibody [MM0014-5F66] (ab51869)

ab51869 at 1/200 dilution, staining PDGF BB in human placental tissue section by Immunohistochemistry (Formalin/PFA fixed paraffin-embedded sections).

Antigen retrieval was done using proteinase K (20 µg/mL in PBS: Room temperature, ranging from 5 min to 30 mins).



Overlay histogram showing SHSY-5Y cells stained with ab51869 (red line). The cells were fixed with 4% paraformaldehyde (10 min) and then permeabilized with 0.1% PBS-Tween for 20 min. The cells were then incubated in 1x PBS / 10% normal goat serum / 0.3M glycine to block non-specific protein-protein interactions followed by the antibody (ab51869, 1 µg/1x10⁶ cells) for 30 min at 22°C. The secondary antibody used was DyLight® 488 goat anti-mouse IgG (H+L) (**ab96879**) at 1/500 dilution for 30 min at 22°C. Isotype control antibody (black line) was mouse IgG1 [ICIGG1] (**ab91353**, 2 µg/1x10⁶ cells) used under the same conditions. Acquisition of >5,000 events was performed. This antibody gave a positive signal in SHSY-5Y cells fixed with 80% methanol (5 min)/permeabilized with 0.1% PBS-Tween for 20 min used under the same conditions.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors