

PE Anti-CD73 antibody [AD2] ab157335

★☆☆☆☆ [1 Abreviews](#) [4 References](#) [1 图像](#)

概述

产品名称	PE Anti-CD73抗体[AD2]
描述	PE小鼠单克隆抗体[AD2] to CD73
宿主	Mouse
偶联物	PE. Ex: 488nm, Em: 575nm
特异性	ab157335 recognizes CD73 expressed predominantly on T and B cell subsets, follicular dendritic cells and endothelial cells
经测试应用	适用于: Flow Cyt
种属反应性	与反应: Human
免疫原	Tissue, cells or virus corresponding to Human CD73. Pre-B leukemia cells
阳性对照	Human blood cells
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 7.4 Preservative: 0.1% Sodium azide Constituents: 0.2% BSA, 99.7% PBS
纯度	Size exclusion
纯化说明	ab157335 is conjugated with R-Phycoerythrin (PE) under optimum conditions. The conjugate is purified by size-exclusion chromatography.
克隆	单克隆
克隆编号	AD2

应用

The Abpromise guarantee

Abpromise™承诺保证使用ab157335于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt	★☆☆☆☆ (1)	Use 10µl for 10 ⁶ cells. Use 10µl per 100µl whole blood. The content of a vial (1 ml) is sufficient for 50 tests. ab91357 - Mouse monoclonal IgG1, is suitable for use as an isotype control with this antibody.

靶标

功能

Hydrolyzes extracellular nucleotides into membrane permeable nucleosides.

疾病相关

Defects in NT5E are the cause of calcification of joints and arteries (CAJA) [MIM:211800]. A condition characterized by adult-onset calcification of the lower extremity arteries, including the iliac, femoral and tibial arteries, and hand and foot capsule joints. Age of onset has been reported as early as the second decade of life, usually involving intense joint pain or calcification in the hands.


序列相似性

Belongs to the 5'-nucleotidase family.

细胞定位

Cell membrane.

图片



Surface staining (flow cytometry) of human peripheral blood with anti-human CD73 PE (ab157335).

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors