

### PE Anti-EVI2B antibody [MEM-216] ab101040

#### 1 图像

#### 概述

产品名称	PE Anti-EVI2B 抗体[MEM-216]
描述	PE 小鼠单克隆抗体[MEM-216] to EVI2B
宿主	Mouse
偶联物	PE. Ex: 488nm, Em: 575nm
经测试应用	适用于: Flow Cyt
种属反应性	与反应: Human
免疫原	Tissue, cells or virus corresponding to Human EVI2B. Raji cells
阳性对照	Flow Cyt: Human peripheral blood cells.
常规说明	<p>The purified antibody is conjugated with R Phycoerythrin (PE) under optimum conditions. The conjugate is purified by size exclusion chromatography and adjusted for direct use. No reconstitution is necessary.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	<p>pH: 7.4</p> <p>Preservative: 0.097% Sodium azide</p> <p>Constituents: 0.2% BSA, PBS</p>
纯度	Protein G purified
纯化说明	>95% as determined by SDS-PAGE
克隆	单克隆
克隆编号	MEM-216

同种型

IgG1

应用

The Abpromise guarantee

Abpromise™承诺保证使用ab101040于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		Use 20µl for 10 <sup>6</sup> cells. <b>ab91357</b> - Mouse monoclonal IgG1, is suitable for use as an isotype control with this antibody.

靶标


组织特异性

Bone marrow, peripheral blood mononuclear cells, fibroblasts and Epstein-Barr virus-transformed lymphoblastoid cell lines.

细胞定位

Membrane.

图片



Surface staining of human peripheral blood using ab101040.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

#### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors