

Anti-Penicillin antibody [P2B9] ab116641

2 References

概述

产品名称	Anti-Penicillin抗体[P2B9]
描述	小鼠单克隆抗体[P2B9] to Penicillin
宿主	Mouse
特异性	ab116641 shows reactivity with Chemical samples.
经测试应用	适用于: IP, ELISA, RIA
种属反应性	与反应: Species independent
免疫原	Penicillin
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at 4°C (up to 6 months). Store at -20°C long term.
存储溶液	Preservative: 0.05% Sodium azide Constituents: 0.02% Potassium chloride, 0.79% Sodium chloride, 0.14% Potassium phosphate, 0.1% BSA
纯度	Protein A purified
克隆	单克隆
克隆编号	P2B9
同种型	IgG1

应用

Abpromise™ 承诺保证使用ab116641于以下的经测试应用

The Abpromise guarantee

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IP		Use at an assay dependent dilution.
ELISA		Use at an assay dependent dilution.
RIA		Use at an assay dependent dilution.

靶标

相关性 Penicillin is a group of Beta-lactam antibiotics used in the treatment of bacterial infections caused by susceptible, usually Gram-positive, organisms. β -lactam antibiotics work by inhibiting the formation of peptidoglycan cross-links in the bacterial cell wall, which results in cytolysis.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors