

PerCP Anti-CD3 antibody [UCHT1], prediluted ab106215

[7 References](#) [1 图像](#)

概述

| | |
|-------|---|
| 产品名称 | PerCP荧光Anti-CD3抗体[UCHT1], prediluted |
| 描述 | PerCP荧光小鼠单克隆抗体[UCHT1] to CD3, prediluted |
| 宿主 | Mouse |
| 偶联物 | PerCP. Ex: 482nm, Em: 675nm |
| 特异性 | ab106215 recognizes the CD3 antigen of the TCR/CD3 complex on mature Human T cells. It reacts with the epsilon chain of the CD3 complex. HLDA I; WS Code T 3 HLDA III; WS Code T 126 HLDA III; WS Code T 471 HLDA VI; WS Code T 6T-CD3.1 |
| 经测试应用 | 适用于: Flow Cyt |
| 种属反应性 | 与反应: Human, Non human primates |
| 免疫原 | Tissue, cells or virus corresponding to Human CD3. Human thymocytes followed by Sezary T cells |
| 阳性对照 | Flow Cyt: Human blood cells. |
| 常规说明 | <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p> |

性能

| | |
|------|---|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Store at +4°C. |
| 存储溶液 | pH: 7.4 Preservative: 0.097% Sodium azide Constituents: 0.2% BSA, 2.68% PBS |
| 纯度 | Size exclusion |
| 克隆 | 单克隆 |
| 克隆编号 | UCHT1 |
| 同种型 | IgG1 |

应用


The Abpromise guarantee

Abpromise™承诺保证使用ab106215于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

| 应用 | Ab评论 | 说明 |
|----------|------|---|
| Flow Cyt | | Use 10µl for 10 ⁶ cells. ab118658 - Mouse monoclonal IgG1, is suitable for use as an isotype control with this antibody. |

图片



Flow Cytometry analysis of Human peripheral blood cells labeling CD3 with Anti-CD3 antibody [UCHT1], prediluted (PerCP) (ab106215).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors