

PerCP/Cy5.5® Anti-Ly-6A/E (Sca-1) antibody [D7] ab93538

1 图像

概述

产品名称	PerCP荧光/Cy5.5® Anti-Ly-6A/E (Sca-1)抗体[D7]
描述	PerCP荧光/Cy5.5®大鼠单克隆抗体[D7] to Ly-6A/E (Sca-1)
宿主	Rat
偶联物	PerCP/Cy5.5®. Ex: 482nm, Em: 690nm
经测试应用	适用于: Flow Cyt
种属反应性	与反应: Mouse
免疫原	Tissue, cells or virus corresponding to Mouse Sca1/ Ly6A/E. IL2 dependent T cell line CTL-L (Mouse).
阳性对照	C57Bl/6 Mouse splenocytes
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	<p>pH: 7.20</p> <p>Preservative: 0.09% Sodium azide</p> <p>Constituents: 0.87% Sodium chloride, 0.1% Gelatin, PBS</p>
纯度	Protein G purified
克隆	单克隆
克隆编号	D7
同种型	IgG2a
轻链类型	kappa

应用

The Abpromise guarantee **Abpromise™** 承诺保证使用 ab93538 于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab 评论	说明
Flow Cyt		

应用说明

Flow Cyt: Use 0.125 - 0.25 µg for 10^5 - 10^8 cells in a 100 µl volume.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

靶标

功能 T-cell activation.

组织特异性 Widely expressed.


序列相似性 Contains 1 UPAR/Ly6 domain.

翻译后修饰 O-glycosylated. Not N-glycosylated.

Not phosphorylated.

细胞定位 Cell membrane.

图片



Flow cytometry staining of C57Bl/6 splenocytes with 0.25 µg of PerCP-Cy5.5 Rat IgG2a Isotype Control (open histogram) or 0.125 µg of ab93538 (colored histogram). Cells in the lymphocyte gate were used for analysis.

Flow Cytometry - PerCP/Cy5.5® Anti-Ly-6A/E (Sca-1) antibody [D7] (ab93538)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors