

Anti-Pneumolysin antibody [PLY-4] ab71810

★★★★★ [1 Abreviews](#) [15 References](#) [1 图像](#)

概述	
产品名称	Anti-Pneumolysin抗体[PLY-4]
描述	小鼠单克隆抗体[PLY-4] to Pneumolysin
宿主	Mouse
经测试应用	适用于: WB
种属反应性	与反应: Species independent
免疫原	Recombinant full length protein
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能	
形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	pH: 7.50 Constituent: 0.06% Tris buffered saline
纯度	Protein G purified
克隆	单克隆
克隆编号	PLY-4
同种型	IgG1

应用	
The Abpromise guarantee	Abpromise™ 承诺保证使用ab71810于以下的经测试应用
“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。	


应用	Ab评论	说明
WB	★★★★★ (1)	Use a concentration of 3 µg/ml.

靶标

细胞定位

Plasma membrane

图片



Coomassie Blue-stained gel (lane 1) and Western blots of purified recombinant PLY, 3 µg per lane, with mouse monoclonal antibodies PLY-4 (lane 2), PLY-5 (lane 3) and PLY-7 (lane 4); a PLY-irrelevant mouse monoclonal antibody 1.4G8.66 served as control (lane 5).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors