

Anti-Salmonella typhi antibody [AC04] ab69238

2 References

概述

产品名称	Anti-Salmonella typhi抗体[AC04]
描述	小鼠单克隆抗体[AC04] to Salmonella typhi
宿主	Mouse
特异性	ab69238 is only reactive with Salmonella typhi.
经测试应用	适用于: ELISA, ICC/IF
种属反应性	与反应: Salmonella typhi
免疫原	Tissue, cells or virus corresponding to Salmonella typhi. Native Salmonella typhi
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: 0.065% Sodium azide Constituent: PBS
纯度	Protein A purified
克隆	单克隆
克隆编号	AC04
骨髓瘤	Sp2/0
同种型	IgG1

应用

The Abpromise guarantee

Abpromise™承诺保证使用ab69238于以下的经测试应用

“应用说明”部分下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		Use at an assay dependent concentration.
ICC/IF		Use at an assay dependent concentration.

靶标

相关性

The genus *Salmonella* is a member of the family Enterobacteriaceae. The genus is composed of Gram-negative bacilli that are facultative and flagellated (motile). Salmonellae possess 3 major antigens; the "H" or flagellar antigen (phase 1 & 2), the "O" or somatic antigen (part of the LPS moiety) and the "Vi" or capsular antigen (referred to as "K" in other Enterobacteriaceae). *Salmonella typhi* is strictly a human pathogen and is the causative agent of Typhoid fever which is characterized by a sustained fever as high as 40°C (104°F), profuse sweating, gastroenteritis, and nonbloody diarrhea

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors