

Anti-Serum Amyloid A antibody [EPR11592-129] ab190801

重组 RabMAb

3 图像

概述

产品名称	Anti-Serum Amyloid A抗体[EPR11592-129]
描述	兔单克隆抗体[EPR11592-129] to Serum Amyloid A
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Synthetic peptide. This information is proprietary to Abcam and/or its suppliers.
阳性对照	Human serum and plasma lysates; Serum Amyloid A recombinant protein (aa1-122 of UniProt P0DJ19, with a proprietary tag).
常规说明	<p>This product is a recombinant monoclonal antibody, which offers several advantages including:</p> <ul style="list-style-type: none"> - High batch-to-batch consistency and reproducibility - Improved sensitivity and specificity - Long-term security of supply - Animal-free production <p>For more information see here.</p> <p>Our RabMAb[®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to RabMAb[®] patents.</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
存储溶液	<p>pH: 7.2</p> <p>Preservative: 0.01% Sodium azide</p> <p>Constituents: 40% Glycerol (glycerin, glycerine), 0.05% BSA, 59% PBS</p>
纯度	Protein A purified
克隆	单克隆
克隆编号	EPR11592-129
同种型	IgG

应用

The Abpromise guarantee

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

Abpromise™承诺保证使用ab190801于以下的经测试应用


应用	Ab评论	说明
WB		1/1000. Detects a band of approximately 12 kDa (predicted molecular weight: 14 kDa).

靶标

相关性

Function: Major acute phase reactant. Apolipoprotein of the HDL complex. Tissue specificity: Expressed by the liver; secreted in plasma. Disease: Note=Reactive, secondary amyloidosis is characterized by the extracellular accumulation in various tissues of the SAA protein. These deposits are highly insoluble and resistant to proteolysis; they disrupt tissue structure and compromise function. Note=Elevated serum SAA protein levels may be associated with lung cancer. Similarity: Belongs to the SAA family. PTM: This protein is the precursor of amyloid protein A, which is formed by the removal of approximately 24 residues from the C-terminal end.

图片



All lanes : Anti-Serum Amyloid A antibody [EPR11592-129] (ab190801) at 1/1000 dilution

Lane 1 : Human serum lysate

Lane 2 : Human plasma lysate


Lysates/proteins at 10 µg per lane.

Secondary

All lanes : Anti-Rabbit IgG (HRP), specific to the non-reduced form of IgG at 1/1000 dilution

Predicted band size: 14 kDa

Observed band size: 12 kDa



Western blot - Anti-Serum Amyloid A antibody
[EPR11592-129] (ab190801)

Anti-Serum Amyloid A antibody [EPR11592-129] (ab190801) at
1/5000 dilution + Human Serum Amyloid A2 recombinant protein
(UniProt P0DJ19) with proprietary tag at 0.01 µg

Secondary

Anti-Rabbit IgG (HRP), specific to the non-reduced form of IgG at
1/1000 dilution

Predicted band size: 14 kDa

Observed band size: 38 kDa

Why choose a recombinant antibody?



Research with confidence
Consistent and reproducible results



Long-term and scalable supply
Recombinant technology



Success from the first experiment
Confirmed specificity



Ethical standards compliant
Animal-free production

Anti-Serum Amyloid A antibody [EPR11592-129]
(ab190801)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors