


Anti-SSEA4 antibody [MC813-70] ab16287

★★★★★ [11 Abreviews](#) [165 References](#) [7 图像](#)

概述

产品名称	Anti-SSEA4抗体[MC813-70]
描述	小鼠单克隆抗体[MC813-70] to SSEA4
宿主	Mouse
经测试应用	适用于: ICC/IF, Flow Cyt
种属反应性	与反应: Human 预测可用于: Mouse, Rat 
免疫原	Tissue, cells or virus corresponding to Human SSEA4. Human embryonal carcinoma cell line 2102Ep

常规说明

This antibody clone is manufactured by Abcam. If you require a custom buffer formulation or conjugation for your experiments, please contact orders@abcam.com.

The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.02% Sodium azide
纯度	Tissue culture supernatant
纯化说明	Tissue culture supernatant was cross flow concentrated.
克隆	单克隆
克隆编号	MC813-70
骨髓瘤	Sp2/0

同种型 IgG3
轻链类型 kappa

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab16287于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。


应用	Ab评论	说明
ICC/IF	★★★★★ (8)	Use at an assay dependent concentration.
Flow Cyt	★☆☆☆☆ (2)	Use at an assay dependent concentration. ab18392 - Mouse monoclonal IgG3, is suitable for use as an isotype control with this antibody.

靶标

相关性 Glycosphingolipids function as mediators of cell adhesion and are modulators of signal transduction. SSEA-4 (Stage-Specific Embryonic Antigen 4) is a glycolipid expressed early in embryonic development and in pluripotent stem cells. This antibody was first described and named as part of a series of embryonic antigens, defined by monoclonal antibodies isolated in the lab of Prof. Davor Solter (Kannagi, R. et al., 1983, EMBO J. 2:2355). SSEA-4 can be used as a marker of human embryonic stem cells, human embryonic carcinoma cells and human embryonic germ cells. Monoclonal antibodies to this target have been widely used in the characterization of pluripotent stem cells. Mouse pluripotent stem cells are not recognised by anti-SSEA-4 antibodies but do express the antigen upon differentiation.

细胞定位 Cell Membrane


图片



Flow Cytometry - Anti-SSEA4 antibody [MC813-70] (ab16287)

This image is courtesy of Professor Chris Denning's lab, University of Nottingham


Human induced pluripotent stem cells (iPSCs) stained with ab16287 (red fill) with secondary only control (blue fill). In brief, iPSCs were fixed in 4% formaldehyde (methanol-free) for 15 min at 25°C. Cells were then incubated with the antibody (ab16287, 1/9600 dilution) for 30 min at 4°C. The secondary antibody used was Alexa Fluor® 488 goat anti-mouse IgG at 1/1000 dilution for 30 min at 4°C. iPSCs differentiated for 14 days toward a cardiomyocyte lineage were used as a negative control. Acquisition of >20,000 total events were collected using a 100mW solid state diode laser (488nm) and 529/28 bandpass filter.



Immunocytochemistry/ Immunofluorescence - Anti-SSEA4 antibody [MC813-70] (ab16287)

This image is courtesy of an Abreview submitted by Sun Jeong Kim


ab16287 staining SSEA4 in rat tendon derived stem cells by ICC/IF (Immunocytochemistry/immunofluorescence). Cells were fixed with paraformaldehyde, permeabilized with 0.25% Triton X-100 in PBS and blocked with 1% BSA for 30 minutes at room temperature. Samples were incubated with primary antibody (1/100) for 20 hours at 4°C. An Alexa Fluor® 594-conjugated donkey anti-mouse IgG polyclonal (1/1000) was used as the secondary antibody.



Immunocytochemistry/ Immunofluorescence - Anti-SSEA4 antibody [MC813-70] (ab16287)

Image from Ferro F et al., PLoS One. 2012;7(7):e41774. Epub 2012 Jul 23. Fig 6.; doi:10.1371/journal.pone.0041774; July 23, 2012, PLoS ONE 7(7): e41774.


Immunofluorescence analysis of Human dental pulp stem cells, staining SSEA4 with ab16287 at 1/40 dilution. A FITC-conjugated anti-mouse IgG was used as the secondary antibody.



Immunocytochemistry/ Immunofluorescence - Anti-SSEA4 antibody [MC813-70] (ab16287)


This image is courtesy of an abreview submitted by Vladimir Milenkovic, University Hospital Regensburg, Germany

Immunocytochemistry/ Immunofluorescence analysis of human iPSC cells labeling SSEA4 with ab16287 at 1/500 dilution. Cells were fixed with paraformaldehyde and permeabilized with 0.5% TX100. Cells were blocked with 5% serum for 20 minutes at 25°C. A goat polyclonal anti-mouse Cy3 secondary antibody at 1/500 dilution was used.




Flow Cytometry - Anti-SSEA4 antibody [MC813-70] (ab16287)

anti-SSEA4 antibody, **ab16297** can be used as a marker of Embryonic Carcinoma cells in Flow Cytometry/FACS. As can be seen from the histograms, in non-differentiating conditions (i.e. without retinoic acid) NTERA2 cells were recognised by **ab16297**. However, upon differentiation (addition of retinoic acid), the antibody lost the ability to recognise the cells.



Immunocytochemistry/ Immunofluorescence - Anti-SSEA4 antibody [MC813-70] (ab16287)

2102ep Human Embryonic Carcinoma cells were stained with SSEA4 antibody ab16287. As expected, staining localised to the cell surface (green). Nuclei are stained blue using Hoechst.



Immunocytochemistry/ Immunofluorescence - Anti-SSEA4 antibody [MC813-70] (ab16287)

This image is a courtesy of Fiona Lewis

ab16287 staining SSEA4 in human Embryonic Stem Cells, HUES7 by Immunocytochemistry/ Immunofluorescence. Cells were fixed with paraformaldehyde, permeabilized with Triton and blocking with 10% serum for 1 hour was performed. Samples were incubated with primary antibody (1/100: in 1% serum, 0.1% Triton in PBS) for 1 hour at 37°C. An Alexa Fluor® 588-conjugated goat polyclonal to mouse IgG was used at dilution at 1/100 as secondary antibody.

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors