

Syrian Hamster IgG [SHG-1] - isotype control ab18426

1 References

概述

产品名称	叙利亚仓鼠IgG [SHG-1] -同型对照
经测试应用	适用于: IHC-Fr, Flow Cyt, IP, WB, IHC-P
免疫原	Trinitrophenol + KLH.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 7.20 Preservative: 0.09% Sodium azide Constituent: PBS
纯度	Protein A purified
纯化说明	This antibody was purified by affinity chromatography.
克隆	单克隆
克隆编号	SHG-1
同种型	IgG

应用

The Abpromise guarantee **Abpromise™** 承诺保证使用ab18426于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-Fr		
Flow Cyt		
IP		
WB		
IHC-P		

应用说明

Flow Cyt: Use at an assay dependent dilution for surface or intracellular staining.
IHC-P: Use at an assay dependent dilution.
IHC-Fr: Use at an assay dependent dilution.
IP: Use at an assay dependent dilution.
WB: Use at an assay dependent dilution.

Not tested in other applications.
Optimal dilutions/concentrations should be determined by the end user.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors