

Anti-Talin 2 antibody [68E7] ab105458

★★★★★ [1 Abreviews](#) [9 References](#)

概述

产品名称	Anti-Talin 2抗体[68E7]
描述	小鼠单克隆抗体[68E7] to Talin 2
宿主	Mouse
经测试应用	适用于: WB, ELISA, ICC/IF, IHC-Fr, IP
种属反应性	与反应: Mouse
免疫原	Recombinant fragment, corresponding to amino acids 489-911 of Mouse Talin 2.
阳性对照	Fibroblasts
常规说明	Growth Conditions: RPMI1650 + Glu + Pen + Strep + 10% FCS @ 37°C + 5%CO2 The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing. If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: 0.02% Sodium azide Constituent: 99.98% PBS
纯度	Protein A purified
克隆	单克隆
克隆编号	68E7
骨髓瘤	NS0/1
同种型	IgG2b

The Abpromise guarantee **Abpromise™**承诺保证使用ab105458于以下的经测试应用

“应用说明”部分下显示的仅为推荐的起始稀释度；实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB	★★★★★ (1)	Use at an assay dependent concentration.
ELISA		Use at an assay dependent concentration.
ICC/IF		Use at an assay dependent concentration.
IHC-Fr		Use at an assay dependent concentration.
IP		Use at an assay dependent concentration.

靶标

功能	As a major component of focal adhesion plaques that links integrin to the actin cytoskeleton, may play an important role in cell adhesion. Recruits PIP5K1C to focal adhesion plaques and strongly activates its kinase activity.
序列相似性	Contains 1 FERM domain. Contains 1 I/LWEQ domain.
细胞定位	Cell junction > focal adhesion. Cell junction > synapse. Cell membrane. Cytoplasm > cytoskeleton. Focal adhesion plaques and synapses.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors