

Anti-TGE Virus peplomer antibody [1Q17] ab2599

概述

产品名称	Anti-TGE Virus peplomer抗体[1Q17]
描述	小鼠单克隆抗体[1Q17] to TGE Virus peplomer
宿主	Mouse
特异性	This antibody is specific for TGE (Transmissible gastroenteritis) virus peplomer antigen.
经测试应用	适用于: WB, ELISA
种属反应性	与反应: Transmissible gastroenteritis
免疫原	Tissue, cells or virus corresponding to TGE Virus peplomer. Purified TGE virus
表位	The exact epitope recognised by this antibody has not been mapped.
常规说明	<p>The peplomer protein mediates the binding of virions to the host cell receptor and is also involved with membrane fusion.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	pH: 7.20 Preservative: 0.01% Sodium azide
Primary antibody说明	The peplomer protein mediates the binding of virions to the host cell receptor and is also involved with membrane fusion.
克隆	单克隆
克隆编号	1Q17
骨髓瘤	unknown

同种型 IgG1
轻链类型 unknown

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab2599于以下的经测试应用

“应用说明”部分下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		1/1000.
ELISA		1/10000.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors