


Anti-TGF beta 2 antibody ab80059

★★★★★ [1 Abreviews](#) [2 References](#)

概述

产品名称	Anti-TGF beta 2抗体
描述	兔多克隆抗体to TGF beta 2
宿主	Rabbit
特异性	Reacts specifically with human 48 kDa TGF beta 2 protein.
经测试应用	适用于: IHC, WB
种属反应性	与反应: Human 预测可用于: a wide range of other species 
免疫原	Synthetic peptide derived from the N terminal domain of Human TGF beta 2.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Constituent: Whole serum
纯度	Whole antiserum
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee **Abpromise™** 承诺保证使用ab80059于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC		1/100 - 1/500.
WB	★★★★★ (1)	1/500 - 1/5000. Predicted molecular weight: 47 kDa.

靶标

功能	TGF-beta 2 has suppressive effects on interleukin-2 dependent T-cell growth.
疾病相关	A chromosomal aberration involving TGFB2 is found in a family with Peters anomaly. Translocation t(1;7)(q41;p21) with HDAC9. Loeys-Dietz syndrome 4 Defects in TGFB2 may be a cause of non-syndromic aortic disease (NSAD). NSAD is a frequently asymptomatic but potentially lethal disease characterized by thoracic aortic aneurysms and dissections without additional syndromic features.
序列相似性	Belongs to the TGF-beta family.
翻译后修饰	The precursor is cleaved into mature TGF-beta-2 and LAP, which remains non-covalently linked to mature TGF-beta-2 rendering it inactive.
细胞定位	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors