

### Anti-TIE2 antibody [Cl. 16] - BSA and Azide free ab24859

★★★★★ [13 Abreviews](#) [24 References](#) [5 图像](#)

#### 概述

产品名称	Anti-TIE2抗体[Cl. 16] - BSA and Azide free
描述	小鼠单克隆抗体[Cl. 16] to TIE2 - BSA and Azide free
宿主	Mouse
经测试应用	适用于: IHC-Fr, WB, Flow Cyt
种属反应性	与反应: Human
免疫原	Recombinant fragment corresponding to Human TIE2. Recombinant human soluble extracellular TIE2. Database link: <a href="#">Q02763</a>
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Constituent: PBS
无载体	是
纯度	Protein G purified
克隆	单克隆
克隆编号	Cl. 16
同种型	IgG1

#### 应用

## The Abpromise guarantee

**Abpromise™** 承诺保证使用ab24859于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-Fr	★★★★★ (2)	Use at an assay dependent concentration.
WB	★★★★★ (1)	Use a concentration of 1 - 2 µg/ml. Predicted molecular weight: 126 kDa.
Flow Cyt		Use 1-2µg for 10 <sup>6</sup> cells. We tested in-house with methanol-fixed cells, but this may not be necessary since the antibody recognises an extracellular epitope. <b>ab170190</b> - Mouse monoclonal IgG1, is suitable for use as an isotype control with this antibody.

## 靶标

### 功能

Tyrosine-protein kinase that acts as cell-surface receptor for ANGPT1, ANGPT2 and ANGPT4 and regulates angiogenesis, endothelial cell survival, proliferation, migration, adhesion and cell spreading, reorganization of the actin cytoskeleton, but also maintenance of vascular quiescence. Has anti-inflammatory effects by preventing the leakage of proinflammatory plasma proteins and leukocytes from blood vessels. Required for normal angiogenesis and heart development during embryogenesis. Required for post-natal hematopoiesis. After birth, activates or inhibits angiogenesis, depending on the context. Inhibits angiogenesis and promotes vascular stability in quiescent vessels, where endothelial cells have tight contacts. In quiescent vessels, ANGPT1 oligomers recruit TEK to cell-cell contacts, forming complexes with TEK molecules from adjoining cells, and this leads to preferential activation of phosphatidylinositol 3-kinase and the AKT1 signaling cascades. In migrating endothelial cells that lack cell-cell adhesions, ANGPT1 recruits TEK to contacts with the extracellular matrix, leading to the formation of focal adhesion complexes, activation of PTK2/FAK and of the downstream kinases MAPK1/ERK2 and MAPK3/ERK1, and ultimately to the stimulation of sprouting angiogenesis. ANGPT1 signaling triggers receptor dimerization and autophosphorylation at specific tyrosine residues that then serve as binding sites for scaffold proteins and effectors. Signaling is modulated by ANGPT2 that has lower affinity for TEK, can promote TEK autophosphorylation in the absence of ANGPT1, but inhibits ANGPT1-mediated signaling by competing for the same binding site. Signaling is also modulated by formation of heterodimers with TIE1, and by proteolytic processing that gives rise to a soluble TEK extracellular domain. The soluble extracellular domain modulates signaling by functioning as decoy receptor for angiopoietins. TEK phosphorylates DOK2, GRB7, GRB14, PIK3R1; SHC1 and TIE1.

### 组织特异性

Detected in umbilical vein endothelial cells. Proteolytic processing gives rise to a soluble extracellular domain that is detected in blood plasma (at protein level). Predominantly expressed in endothelial cells and their progenitors, the angioblasts. Has been directly found in placenta and lung, with a lower level in umbilical vein endothelial cells, brain and kidney.

### 疾病相关

Dominantly inherited venous malformations  
May play a role in a range of diseases with a vascular component, including neovascularization of tumors, psoriasis and inflammation.

### 序列相似性

Belongs to the protein kinase superfamily. Tyr protein kinase family. Tie subfamily.  
Contains 3 EGF-like domains.

Contains 3 fibronectin type-III domains.  
Contains 2 Ig-like C2-type (immunoglobulin-like) domains.  
Contains 1 protein kinase domain.

## 结构域

The soluble extracellular domain is functionally active in angiopoietin binding and can modulate the activity of the membrane-bound form by competing for angiopoietins.

## 翻译后修饰


Proteolytic processing leads to the shedding of the extracellular domain (soluble TIE-2 alias sTIE-2).


Autophosphorylated on tyrosine residues in response to ligand binding. Autophosphorylation occurs in trans, i.e. one subunit of the dimeric receptor phosphorylates tyrosine residues on the other subunit. Autophosphorylation occurs in a sequential manner, where Tyr-992 in the kinase activation loop is phosphorylated first, followed by autophosphorylation at Tyr-1108 and at additional tyrosine residues. ANGPT1-induced phosphorylation is impaired during hypoxia, due to increased expression of ANGPT2. Phosphorylation is important for interaction with GRB14, PIK3R1 and PTPN11. Phosphorylation at Tyr-1102 is important for interaction with SHC1, GRB2 and GRB7. Phosphorylation at Tyr-1108 is important for interaction with DOK2 and for coupling to downstream signal transduction pathways in endothelial cells. Dephosphorylated by PTPRB. Ubiquitinated. The phosphorylated receptor is ubiquitinated and internalized, leading to its degradation.

## 细胞定位

Cell membrane. Cell junction. Cell junction, focal adhesion. Cytoplasm, cytoskeleton. Secreted. Recruited to cell-cell contacts in quiescent endothelial cells. Colocalizes with the actin cytoskeleton and at actin stress fibers during cell spreading. Recruited to the lower surface of migrating cells, especially the rear end of the cell. Proteolytic processing gives rise to a soluble extracellular domain that is secreted.


## 图片





Immunohistochemistry (Frozen sections) - Anti-TIE2 antibody [Cl. 16] - BSA and Azide free (ab24859)

ab24859 staining TIE2 in Human spleen by Immunohistochemistry (Frozen sections).




Immunohistochemistry (Frozen sections) - Anti-TIE2 antibody [Cl. 16] - BSA and Azide free (ab24859)

This image is courtesy of an anonymous Abreview

Immunohistochemical analysis of Human brain tissue, staining TIE2 with ab24859.


Tissue was fixed with paraformaldehyde, permeabilized with 0.25 Triton X-100 and blocked with 2.5% BSA for 30 minutes at 25°C. Samples were incubated with primary antibody (1/200 in 2.5% horse serum) for 18 hours at 4°C. An HRP-conjugated horse anti-rabbit polyclonal IgG was used as the secondary antibody.



Flow Cytometry - Anti-TIE2 antibody [Cl. 16] - BSA and Azide free (ab24859)

Overlay histogram showing JEG-3 cells stained with ab24859 (red line). The cells were fixed with methanol (5 min) and incubated in 1x PBS / 10% normal goat serum / 0.3M glycine to block non-specific protein-protein interactions. The cells were then incubated with the antibody (ab24859, 2µg/1x10<sup>6</sup> cells) for 30 min at 22°C. The secondary antibody used was DyLight® 488 goat anti-mouse IgG (H+L) ([ab96879](#)) at 1/500 dilution for 30 min at 22°C. Isotype control antibody (black line) was mouse IgG1 [ICIGG1] ([ab91353](#), 2µg/1x10<sup>6</sup> cells) used under the same conditions. Acquisition of >5,000 events was performed. This antibody gave a decreased signal in JEG-3 cells fixed with 4% paraformaldehyde (10 min) used under the same conditions.

Please note that Abcam does not have data for use of this antibody on non-fixed cells. We welcome any customer feedback.



**All lanes :** Anti-TIE2 antibody [Cl. 16] - BSA and Azide free (ab24859)

**Lane 1 :** HUVECs left untreated

**Lane 2 :** HUVECs stimulated for 3 hours with PMA at 25 ng/ml

**Lane 3 :** HUVECs stimulated for 6 hours with PMA at 25 ng/ml

**Lane 4 :** HUVECs stimulated for 9 hours with PMA at 25 ng/ml

**Lane 5 :** HUVECs stimulated for 24 hours with PMA at 25 ng/ml

**Predicted band size:** 126 kDa

Samples were immunoprecipitated with another TIE2 antibody.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

#### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors