

### Anti-TLR3 antibody ab62567

#### 1 图像

#### 概述

|       |   |
|-------|---|
| 产品名称  | Anti-TLR3抗体   |
| 描述    | 兔多克隆抗体to TLR3   |
| 宿主    | Rabbit  |
| 经测试应用 | 适用于: WB   |
| 种属反应性 | 与反应: Mouse  |
| 免疫原   | Synthetic peptide within Human TLR3 aa 40-90 (N terminal). The exact sequence is proprietary.<br>15 aa long.<br>Database link: <a href="#">O15455</a>   |
| 阳性对照  | WB: EL4 cell lysate.  |
| 常规说明  | <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p> |

#### 性能

|      |   |
|------|---|
| 形式   | Liquid  |
| 存放说明 | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle. |
| 存储溶液 | pH: 7.2<br>Preservative: 0.02% Sodium azide<br>Constituent: PBS   |
| 纯度   | Ion Exchange Chromatography   |
| 克隆   | 多克隆   |
| 同种型  | IgG   |

#### 应用

The Abpromise guarantee


Abpromise™承诺保证使用ab62567于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

| 应用 | Ab评论 | 说明  |
|----|------|---|
| WB |      | Use a concentration of 1 - 2 µg/ml. Detects a band of approximately 90 kDa (predicted molecular weight: 104 kDa). |

|       |   |
|-------|---|
| 靶标    |   |
| 功能    | Key component of innate and adaptive immunity. TLRs (Toll-like receptors) control host immune response against pathogens through recognition of molecular patterns specific of microorganisms. TLR3 is a nucleotide-sensing TLR which is activated by double-stranded RNA, a sign of viral infection. Acts via MYD88 and TRAF6, leading to NF-kappa-B activation, cytokine secretion and the inflammatory response.   |
| 组织特异性 | Expressed at high level in placenta and pancreas. Also detected in CD11c+ immature dendritic cells. Only expressed in dendritic cells and not in other leukocytes, including monocyte precursors. TLR3 is the TLR that is expressed most strongly in the brain, especially in astrocytes, glia, and neurons.  |
| 疾病相关  | Defects in TLR3 are associated with herpes simplex encephalitis type 2 (HSE2) [MIM:613002]. HSE is a rare complication of human herpesvirus 1 (HHV-1) infection, occurring in only a small minority of HHV-1 infected individuals. HSE is characterized by hemorrhagic necrosis of parts of the temporal and frontal lobes. Onset is over several days and involves fever, headache, seizures, stupor, and often coma, frequently with a fatal outcome. Note=TLR3 mutations predispose otherwise healthy individuals to isolated herpes simplex encephalitis through a mechanism that involves impaired IFNs production and reduced immune defense against viral infection in the central nervous system. |
| 序列相似性 | Belongs to the Toll-like receptor family.<br>Contains 22 LRR (leucine-rich) repeats.<br>Contains 1 LRRCT domain.<br>Contains 1 LRRNT domain.<br>Contains 1 TIR domain.  |
| 结构域   | ds-RNA binding is mediated by LRR 1 to 3, and LRR 17 to 18.   |
| 翻译后修饰 | Heavily N-glycosylated, except on that part of the surface of the ectodomain that is involved in ligand binding.  |
| 细胞定位  | Endoplasmic reticulum membrane. Endosome membrane.  |

图片



Anti-TLR3 antibody (ab62567) at 0.5 µg/ml + EL4 (mouse lymphoma cell line) cell lysate

**Predicted band size:** 104 kDa

**Observed band size:** 90 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

#### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors