


Anti-TMLHE antibody ab97842

1 图像

概述

产品名称	Anti-TMLHE抗体
描述	兔多克隆抗体to TMLHE
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Human 预测可用于: Mouse, Rat, Cow 
免疫原	Recombinant fragment corresponding to Human TMLHE aa 100-400. NP_060666 Database link: Q9NVH6-1
阳性对照	HeLaS3 whole cells lysate; HepG2 and H1299 cells
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	pH: 7.00 Preservative: 0.01% Thimerosal (merthiolate) Constituents: 1.21% Tris, 0.75% Glycine, 10% Glycerol (glycerin, glycerine)
纯度	Protein A purified
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee


Abpromise™承诺保证使用ab97842于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		1/100 - 1/1000. Predicted molecular weight: 50 kDa.

靶标	
功能	Converts trimethyllysine (TML) into hydroxytrimethyllysine (HTML).
组织特异性	All isoforms, but isoform 8, are widely expressed in adult and fetal tissues. Isoform 8 is restricted to heart and skeletal muscle.
通路	Amine and polyamine biosynthesis; carnitine biosynthesis.
疾病相关	Autism, X-linked 6
序列相似性	Belongs to the gamma-BBH/TMLD family.
细胞定位	Mitochondrion matrix.

图片



Western blot - Anti-TMLHE antibody (ab97842)

Anti-TMLHE antibody (ab97842) at 1/200 dilution + HeLaS3 whole cell lysate at 30 µg

Predicted band size: 50 kDa

12% SDS PAGE

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors