


Anti-TPN antibody ab13518

★★★★★ [1 Abreviews](#) [7 References](#)

概述

产品名称	Anti-TPN抗体
描述	兔多克隆抗体to TPN
宿主	Rabbit
特异性	Weakly cross-reacts with bovine Tapasin.
经测试应用	适用于: WB, IP
种属反应性	与反应: Human 预测可用于: Cow 
免疫原	Synthetic peptide corresponding to Human TPN aa 430-448. Sequence: LGWAAVYLSTCKDSKKKAE

 [Run BLAST with](#)

 [Run BLAST with](#)

常规说明

The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.09% Sodium azide Constituents: PBS, 50% Glycerol (glycerin, glycerine)
纯度	Protein A purified
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee

Abpromise™承诺保证使用ab13518于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB	★★★★★ (1)	Use at an assay dependent concentration. Predicted molecular weight: 48 kDa.
IP		Use at an assay dependent concentration.

靶标

功能	Involved in the association of MHC class I with transporter associated with antigen processing (TAP) and in the assembly of MHC class I with peptide (peptide loading).
组织特异性	Neutrophils, mostly in fully differentiated cells.
序列相似性	Contains 1 Ig-like C1-type (immunoglobulin-like) domain.
结构域	The N-terminus is required for efficient association with MHC class I molecule and for a stable interaction between MHC I and calreticulin. Binding to TAP is mediated by the C-terminus region.
细胞定位	Endoplasmic reticulum membrane.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors