

# Anti-Urease antibody ab34806

## 1 References

### 概述

|       |   |
|-------|---|
| 产品名称  | Anti-Urease抗体   |
| 描述    | 兔多克隆抗体to Urease   |
| 宿主    | Rabbit  |
| 经测试应用 | 适用于: ELISA  |
| 种属反应性 | 与反应: Canavalia ensiformis   |
| 免疫原   | Full length protein corresponding to Urease. Urease (Jack Bean)   |
| 常规说明  | <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p> |

### 性能

|      |  |
|------|--|
| 形式   | Liquid   |
| 存放说明 | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.            |
| 存储溶液 | Preservative: 0.01% Sodium azide<br>Constituents: 0.88% Sodium chloride, 0.164% Sodium phosphate |
| 纯度   | Whole antiserum  |
| 克隆   | 多克隆  |
| 同种型  | IgG  |

### 应用

#### The Abpromise guarantee

**Abpromise™**承诺保证使用ab34806于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

| 应用    | Ab评论 | 说明                                       |
|-------|------|--|
| ELISA |      | Use at an assay dependent concentration. |

## 靶标

|      |   |
|------|---|
| 相关性  | The enzyme urease catalyzes the hydrolysis of urea to create carbon dioxide and ammonia. The protein consists of a single kind of polypeptide chain containing 840 amino acid residues. The subunit relative molecular mass calculated from the sequence is 90,770, indicating that urease is composed of six subunits. |
| 细胞定位 | Cytoplasmic   |

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors