

Anti-Uteroglobin antibody ab40873

★★★★★ [6 Abreviews](#) [37 References](#) [2 图像](#)

概述

产品名称	Anti-Uteroglobin抗体
描述	兔多克隆抗体to Uteroglobin
宿主	Rabbit
经测试应用	适用于: IHC-P
种属反应性	与反应: Mouse
免疫原	Recombinant rat uteroglobin His fusion protein.
常规说明	

This antibody is useful for distinguishing ciliated from non-ciliated cells in the adult mouse lung.

The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: 0.01% Sodium azide Constituent: Whole serum
纯度	Whole antiserum
Primary antibody说明	This antibody is useful for distinguishing ciliated from non-ciliated cells in the adult mouse lung.
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用 ab40873 于以下的经测试应用

“应用说明”部分下显示的仅为推荐的起始稀释度；实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-P	★★★★★ (3)	1/1600.

靶标


功能	Binds phosphatidylcholine, phosphatidylinositol, polychlorinated biphenyls (PCB) and weakly progesterone, potent inhibitor of phospholipase A2.
组织特异性	Clara cells (nonciliated cells of the surface epithelium of the pulmonary airways).
序列相似性	Belongs to the secretoglobin family.
细胞定位	Secreted.

图片



Adult mouse lung was embedded in paraffin and stained with ab40873 (1:1600 dilution).

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Uteroglobin antibody (ab40873)



ab40873 staining Uteroglobin in human lung parenchyma tissue by Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections). Tissue was fixed with formaldehyde and a heat mediated antigen retrieval step was performed using 10mM citrate buffer. Samples were then blocked with 1% BSA, 0.1% Triton-X/ PBS for 15 minutes at 18°C followed by incubation with the primary antibody at a 1/500 dilution for 16 hours at 4°C. An HRP-conjugated pig anti-rabbit polyclonal was used as secondary antibody at a 1/100 dilution.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Uteroglobin antibody (ab40873)

Image courtesy of an anonymous Abreview.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors