




Anti-V5 tag antibody [SV5-P-K] ab206565

重组

2 图像

概述

产品名称	Anti-V5 tag抗体[SV5-P-K]
描述	仓鼠单克隆抗体[SV5-P-K] to V5 tag
宿主	Hamster
特异性	ab206565 recognises a 14 amino acid sequence (GKPIPPLLGLDST) derived from the P and V proteins of the Paramyxo Virus simian virus 5.
经测试应用	适用于: WB, ELISA, IHC-Fr, ICC/IF, IP, RIA
种属反应性	与反应: Species independent
免疫原	Synthetic peptide corresponding to V5 tag. Paramyxo Virus Simian Virus 5.  <a href="#">Run BLAST with</a>  <a href="#">Run BLAST with</a> 

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.02% Proclin 300 Constituent: 99% PBS
纯度	Protein A purified
克隆	单克隆
克隆编号	SV5-P-K
同种型	IgG
轻链类型	kappa

应用

The Abpromise guarantee      **Abpromise™**承诺保证使用ab206565于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。


应用	Ab评论	说明
WB		Use at an assay dependent concentration.
ELISA		Use at an assay dependent concentration.
IHC-Fr		Use at an assay dependent concentration.
ICC/IF		Use at an assay dependent concentration.
IP		Use at an assay dependent concentration.
RIA		Use at an assay dependent concentration.

## 靶标

### 相关性

The V5 epitope tag is derived from a small epitope (Pk) present on the P and V proteins of the paramyxovirus of simian virus 5 (SV5). The V5 tag is usually used with all 14 amino acids (GKPIPNPLLGLDST), although it has also been used with a shorter 9 amino acid sequence (IPNPLLGLD).

## 图片



ELISA - Anti-V5 tag antibody [SV5-P-K] (ab206565)

Plate coated with 5 µg/ml of multi-tag protein. ab206565 loaded at 1000 ng/ml in first well and an 11 point serial dilution performed in duplicate for each antibody with last well left blank. Signal detection was performed using an HRP-conjugated goat anti-mouse.

### Why choose a recombinant antibody?



**Research with confidence**  
Consistent and reproducible results



**Long-term and scalable supply**  
Recombinant technology



**Success from the first experiment**  
Confirmed specificity



**Ethical standards compliant**  
Animal-free production

Anti-V5 tag antibody [SV5-P-K] (ab206565)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors