

Anti-V5 tag antibody [SV5-P-K] ab206571

重组

2 References 2 图像

概述

产品名称	Anti-V5 tag抗体[SV5-P-K]
描述	大鼠单克隆抗体[SV5-P-K] to V5 tag
宿主	Rat
特异性	ab206571 recognises a 14 amino acid sequence (GKPIPPLLGLDST) derived from the P and V proteins of the Paramyxo Virus simian virus 5.
经测试应用	适用于: ICC/IF, IP, RIA, WB, ELISA, IHC-Fr
种属反应性	与反应: Species independent
免疫原	Synthetic peptide corresponding to V5 tag aa 50-150. Paramyxo Virus Simian Virus 5. Database link: P11207
常规说明	This is a chimeric antibody created for improved compatibility with existing reagents, assays and techniques.

 [Run BLAST with](#)

 [Run BLAST with](#)

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.02% Proclin 300 Constituent: 99% PBS
纯度	Protein A purified
克隆	单克隆
克隆编号	SV5-P-K
同种型	IgG2a
轻链类型	kappa

应用

Abpromise™ 承诺保证使用ab206571于以下的经测试应用

The Abpromise guarantee

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC/IF		Use at an assay dependent concentration.
IP		Use at an assay dependent concentration.
RIA		Use at an assay dependent concentration.
WB		Use at an assay dependent concentration.
ELISA		Use at an assay dependent concentration.
IHC-Fr		Use at an assay dependent concentration.

靶标

相关性

The V5 epitope tag is derived from a small epitope (Pk) present on the P and V proteins of the paramyxovirus of simian virus 5 (SV5). The V5 tag is usually used with all 14 amino acids (GKPIPPLLGLDST), although it has also been used with a shorter 9 amino acid sequence (IPNPLLGLD).

图片





Plate coated with 5ug/ml of multi-tag protein.

ab206751 loaded at 1000 ng/ml in first well and an 11 point serial dilution performed in duplicate for each antibody with last well left blank. Signal detection was performed using HRP-conjugated goat anti-mouse.

Why choose a recombinant antibody?



Research with confidence
Consistent and reproducible results

Long-term and scalable supply
Recombinant technology

Success from the first experiment
Confirmed specificity

Ethical standards compliant
Animal-free production

Anti-V5 tag antibody [SV5-P-K] (ab206571)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors