

### ACE2 peptide ab198988

[1 References](#) [1 图像](#)

#### 描述

产品名称	ACE2多肽
Accession	<a href="#">Q9BYF1</a>
无动物成分	No
性质	Synthetic

#### 技术指标

Our [Abpromise guarantee](#) covers the use of **ab198988** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Blocking - Blocking peptide for Anti-ACE2 antibody [EPR4435(2)] ( <a href="#">ab108252</a> )
形式	Lyophilized
补充说明	<p>This is the blocking peptide for <a href="#">ab108252</a></p> <ul style="list-style-type: none"><li>- First try to dissolve a small amount of peptide in either water or buffer. The more charged residues on a peptide, the more soluble it is in aqueous solutions.</li><li>- If the peptide doesn't dissolve try an organic solvent e.g. DMSO, then dilute using water or buffer.</li><li>- Consider that any solvent used must be compatible with your assay. If a peptide does not dissolve and you need to recover it, lyophilise to remove the solvent.</li><li>- Gentle warming and sonication can effectively aid peptide solubilisation. If the solution is cloudy or has gelled the peptide may be in suspension rather than solubilised.</li><li>- Peptides containing cysteine are easily oxidised, so should be prepared in solution just prior to use.</li></ul>

#### 制备和贮存

稳定性和存储	Shipped at 4°C. Store at -20°C. Information available upon request.
--------	--

#### 常规信息

功能	Carboxypeptidase which converts angiotensin I to angiotensin 1-9, a peptide of unknown function,
----	--

and angiotensin II to angiotensin 1-7, a vasodilator. Also able to hydrolyze apelin-13 and dynorphin-13 with high efficiency. May be an important regulator of heart function. In case of human coronaviruses SARS and HCoV-NL63 infections, serve as functional receptor for the spike glycoprotein of both coronaviruses.

#### 组织特异性

Expressed in endothelial cells from small and large arteries, and in arterial smooth muscle cells. Expressed in lung alveolar epithelial cells, enterocytes of the small intestine, Leydig cells and Sertoli cells (at protein level). Expressed in heart, kidney, testis, and gastrointestinal system.

#### 序列相似性

Belongs to the peptidase M2 family.

#### 翻译后修饰

N-glycosylation on Asn-90 may limit SARS infectivity.

#### 细胞定位

Secreted and Cell membrane.

## 图片

**Choose the best for your experiment from thousands of proteins and peptides**

- Reinforce the validity of your results with **protein controls**
- Validate specific, reliable reagents with **blocking peptides**
- Keep working on targets without specific antibodies with **fusion-tagged proteins**
- Get consistent, reproducible results with premium **bioactive proteins**

ACE2 peptide (ab198988)

To learn more about our protein and peptide range click [here](#).

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors