

# Human coronavirus SARS-CoV-2 spike glycoprotein peptide ab273063

### 3 References

#### 描述

---

产品名称	人冠状病毒SARS-CoV-2 spike glycoprotein多肽
Accession	<b><u>QHD43416</u></b>
无动物成分	No
性质	Synthetic
种属	Human coronavirus
额外的序列信息	(C terminal). 20 amino acid peptide within the last 50 amino acids of the SARS CoV-2 spike glycoprotein

#### 技术指标

---

Our **Abpromise guarantee** covers the use of **ab273063** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Blocking - Blocking peptide for Anti-SARS-CoV-2 spike glycoprotein antibody - Coronavirus ( <b><u>ab272504</u></b> )
形式	Liquid

#### 制备和贮存

---

稳定性和存储	Shipped at 4°C. Store at -20°C. pH: 7.20 Preservative: 0.02% Sodium azide Constituents: PBS, 0.1% BSA
--------	--

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors