

Human proCNP peptide ab42244

描述	
产品名称	人proCNP多肽
纯度	> 95 % HPLC.
无动物成分	No
性质	Synthetic
种属	Human
序列	GDKAPGGGGANLKGDRSRLLR
氨基酸	30 to 50
技术指标	
Our Abpromise guarantee covers the use of ab42244 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	Blocking
形式	Liquid
制备和贮存	
稳定性和存储	Shipped at 4°C. After reconstitution store at -20°C. Avoid freeze / thaw cycles. Constituent: 0.1% Acetic acid
常规信息	
功能	Hormone which plays a role in endochondral ossification through regulation of cartilaginous growth plate chondrocytes proliferation and differentiation. May also be vasoactive and natriuretic. Specifically binds and stimulates the cGMP production of the NPR2 receptor. Binds the clearance receptor NPR3.
序列相似性	Belongs to the natriuretic peptide family.
细胞定位	Secreted.

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors