

Recombinant Cytomegalovirus ICP36 protein ab43044

描述

产品名称	重组Cytomegalovirus ICP36蛋白
生物活性	Immunoreactive with sera of CMV infected individuals, excellent antigen for detection of CMV with minimal specificity problems.
纯度	> 95 % SDS-PAGE. Sephadex derived purification.
表达系统	Escherichia coli
蛋白长度	Protein fragment
无动物成分	No
性质	Recombinant
额外的序列信息	The protein contains the CMV ICP36 immunodominant regions.

技术指标

Our **Abpromise guarantee** covers the use of **ab43044** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	ELISA Western blot
形式	Liquid

制备和贮存

稳定性和存储	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. Constituents: 0.307% Glutathione, 0.79% Tris HCl, 50% Glycerol, 0.348% Sodium chloride
--------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

常规信息

相关性	In extracts from CMV-infected cells, the CMV DNA polymerase is found strongly associated with the additional polypeptide, ICP36. This protein has been identified as the CMV homolog of the herpes simplex virus type 1 UL42 gene product and may have a similar function. HCMV DNA polymerase and ICP36 expressed in the same system interact to form a stable complex. Moreover, ICP36 functions to stimulate the DNA polymerase activity in a template-dependent
-----	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

manner.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors