

# Recombinant Cytomegalovirus pp28 protein ab43038

## 描述

产品名称	重组Cytomegalovirus pp28蛋白
生物活性	Immunoreactive with sera of CMV infected individuals.
纯度	> 95 % SDS-PAGE. Sephadex derived purification.
表达系统	Escherichia coli
蛋白长度	Protein fragment
无动物成分	No
性质	Recombinant
种属	Cytomegalovirus
额外的序列信息	The protein contains the CMV Pp28 (UL99) immunodominant regions.

## 技术指标

Our **Abpromise guarantee** covers the use of **ab43038** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	ELISA Western blot
形式	Liquid

## 制备和贮存

稳定性和存储	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. pH: 7.20 Constituents: 0.395% Tris HCl, 0.0292% EDTA, 50% Glycerol
--------	--

## 常规信息

相关性	Cytomegalovirus (CMV), is a genus of Herpes viruses; in humans the species is known as Human herpesvirus 5 (HHV-5). It belongs to the Betaherpesvirinae subfamily of Herpesviridae. The name means "very big cell virus".
-----	---



**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors