

Recombinant Dog Anionic Trypsin protein (His tag) ab225607

描述

产品名称	重组狗Anionic Trypsin蛋白(His tag)
纯度	> 90 % SDS-PAGE.
表达系统	Escherichia coli
Accession	<u>P06872</u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Dog
序列	IVGGYTCEENSVPYQVSLNAGYHFCGGLISDQWVSAAHCY KSRIQVRL GEYNIDVLEGNEQFINSKAVIRHPNYNSWILDNDIMLIKLS PAVLNARV ATISLPRACAAPGTQCLISGWGNTLSSGTNYPELLQCLDAPI LTQAQCEA SYPGQITENMICAGFLEGGKDSCQGDSSGGPVVVCNGELQGIVS WYGCAQK NKPGVYTKVCNFVDWIQSTIAANS
预测分子量	40 kDa including tags
氨基酸	24 to 247
标签	His tag N-Terminus
额外的序列信息	This product is the full-length mature protein from aa 24 to 247 with a 6xHis-SUMO tag at the N-terminus. The signal peptide and the propeptide are not included.

技术指标

Our **Abpromise guarantee** covers the use of **ab225607** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE
形式	Liquid

制备和贮存

稳定性和存储

Shipped at 4°C. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 7.2

Constituents: 50% Glycerol (glycerin, glycerine), Tris buffer

常规信息

序列相似性

Belongs to the peptidase S1 family.

Contains 1 peptidase S1 domain.

细胞定位

Secreted, extracellular space.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors