

Recombinant dog Urokinase protein ab92592

1 References 1 图像

描述	
产品名称	重组狗Urokinase蛋白
纯度	> 95 % SDS-PAGE. ab92592 is > 98% pure by SDS-PAGE.
表达系统	Insect cells
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Dog
序列	SHEFHPVSSASNCSLNGGTCVSYKYFSNIQRCNCPKKFQGE HCEIDTSK TCYEGNGHSYRGKANTNTRGRPCLAWNSAAVLLKQYHALRPD ALQLGLGK HNYCRNPDHQRKPWCYVQVGLMQLVQECTVQHCSFEKNPHSP PEKAEFQC GQKALRPRFKIIGGEFTTIENQPWFAAIYRKHRGGSVTYQCG GSLISPCW
氨基酸	1 to 543

技术指标	
Our <b>Abpromise guarantee</b> covers the use of <b>ab92592</b> in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	Functional Studies SDS-PAGE
形式	Liquid
补充说明	Produced using non-baculovirus insect cells.

制备和贮存	
稳定性和存储	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 5.00


Constituents: 0.41% Sodium acetate, 0.0292% EDTA, 0.58% Sodium chloride

This product is an active protein and may elicit a biological response in vivo, handle with caution.

## 常规信息

功能	Specifically cleave the zymogen plasminogen to form the active enzyme plasmin.
组织特异性	Expressed in the prostate gland and prostate cancers.
序列相似性	Belongs to the peptidase S1 family. Contains 1 EGF-like domain. Contains 1 kringle domain. Contains 1 peptidase S1 domain.
翻译后修饰	Phosphorylation of Ser-158 and Ser-323 abolishes proadhesive ability but does not interfere with receptor binding.
细胞定位	Secreted.

## 图片



### 10% SDS-PAGE

Lane 1: Urokinase (3 ug) Reduced

Lane 2: Urokinase (3 ug) + CPAI (12 ug) Reduced

Lane 3: Urokinase (3 ug) Non-reduced

Lane 4: Urokinase (3 ug) + CPAI (12 ug) Non-reduced

Lane 5: Prestained Standard

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors

