

Recombinant Hepatitis E Virus ORF2 antigen protein ab49036

描述	
产品名称	重组Hepatitis E Virus ORF2 antigen蛋白
纯度	> 95 % SDS-PAGE.
表达系统	Escherichia coli
蛋白长度	Protein fragment
无动物成分	No
性质	Recombinant
序列	SANGEPTVKLYTSVENAQQDKGIAIPHDIDLGESRVVIQDYD NQHEQDRP TPSPAPSRP
氨基酸	403 to 461
技术指标	
Our Abpromise guarantee covers the use of ab49036 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	Western blot ELISA
形式	Liquid
制备和贮存	
稳定性和存储	Shipped at 4°C. Store at -20°C. Avoid freeze / thaw cycle. pH: 7.2 Constituents: 48% Urea, 0.078% Beta mercaptoethanol, 0.316% Tris HCl
常规信息	
相关性	The Hepatitis E Virus ORF2 antigen contains a high basic amino acid content, suggesting that it may be involved in the encapsidation of the genomic RNA by effectively neutralizing the negatively charged RNA.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors