

## Product datasheet

# Recombinant HPV16 E6 protein (His tag) ab226447

[2 References](#) [3 图像](#)

### 描述

产品名称	重组HPV16 E6蛋白(His tag)	
纯度	> 90 % SDS-PAGE.	
表达系统	Escherichia coli	
Accession	<b><u>P03126</u></b>	
蛋白长度	Full length protein	
无动物成分	No	
性质	Recombinant	
序列	MHQKRTAMFQDPQERPRKLPQLCTELQTTIHDIILECVYCKQ QLLRREVY DFAFRDLCIVYRDGNPYAVCDKCLKFYISKISEYRHYCYSLYG TTLEQQYN KPLCDLLIRCINCQKPLCPEEKQRHLDKKQRFHNIIRGRWTGR CMSCCRSS RTRRETQL	
预测分子量	23 kDa including tags	
氨基酸	1 to 158	
标签	His tag N-Terminus	
额外的序列信息	Human papillomavirus type 16	

### 技术指标

Our **Abpromise guarantee** covers the use of **ab226447** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE
	Mass Spectrometry
质谱法	LC-MS/MS
形式	Liquid


### 制备和贮存

**稳定性和存储** Shipped at 4°C. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 7.20


Constituents: Tris buffer, 50% Glycerol (glycerin, glycerine)

## 图片




SDS-PAGE - Recombinant HPV16 E6 protein (His tag) (ab226447)

(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Mass Spectrometry - Recombinant HPV16 E6 protein (His tag) (ab226447)

Based on the SEQUEST from database of E.coli host and target protein, the LC-MS/MS Analysis result of ab226447 could indicate that this peptide derived from E.coli-expressed Homo sapiens (Human) E6.



Mass Spectrometry - Recombinant HPV16 E6 protein (His tag) (ab226447)

Based on the SEQUEST from database of E.coli host and target protein, the LC-MS/MS Analysis result of ab226447 could indicate that this peptide derived from E.coli-expressed Homo sapiens (Human) E6.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors