

Recombinant Human ADAMDEC1 protein ab191930

描述

产品名称	重组人ADAMDEC1蛋白
纯度	> 95 % SDS-PAGE. Purity is greater than 95% as determined by SEC-HPLC and reducing SDS-PAGE.
内毒素水平	< 1.000 Eu/μg
表达系统	HEK 293 cells
Accession	<u>O15204</u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	KSPEKEDFLRAQKYIDLYLVLDNAFYKKNYNENLTLIRSFVFD VMNLLNVI YNTIDVQVALVGMEIWSGDGKIKVVPASSTTFDNFLRWSSN LGKKIHDH AQLLSGISFNRRRVGLAASNLCSPSSVAVIEAKKKNVALV GVMSHELG HVLGMPDVPFNTKCPGSCVMNQYLSSKFPKDFSTSCRAHFE RYLLSQKP KLLQAPIPTNIMTTPVCGNHLLVGEDCDCGSPKECTNLCC EALTCKLK PGTDCGGDAPNHTTEVDHHHHHH
预测分子量	51 kDa including tags
氨基酸	206 to 470
标签	His tag C-Terminus

技术指标

Our **Abpromise guarantee** covers the use of **ab191930** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE HPLC
形式	Liquid

制备和贮存

稳定性和存储

Shipped on Dry Ice. The lyophilized protein is stable for a few weeks at room temperature. Store at -20°C or -80°C. Avoid freeze / thaw cycle. Reconstitute for long term storage.

pH: 7.50

Constituents: 0.32% Tris HCl, 0.87% Sodium chloride, 0.01% Zinc chloride

常规信息

功能

May play an important role in the control of the immune response and during pregnancy.

组织特异性

Expressed highly in the small intestine and appendix, moderately in lymph node, mucosal lining of the colon, thymus, spleen and very weakly in the bone marrow. Predominantly expressed in dendritic cells (DC) of the germinal center. Weakly expressed in monocyte and highly expressed in macrophage. Absent in immature DC.

序列相似性

Contains 1 disintegrin domain.

Contains 1 peptidase M12B domain.

细胞定位

Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors