

Recombinant human CD30 protein ab157290

描述

产品名称	重组人CD30蛋白
生物活性	Binds Human CD30L (CD153).
纯度	> 95 % SDS-PAGE.

内毒素水平	< 0.100 Eu/μg
表达系统	HEK 293 cells
Accession	<u>P28908</u>
蛋白长度	Protein fragment
无动物成分	No
性质	Recombinant
种属	Human
序列	

```

MRVLLAALGL LFLGALRAFP QDRPFEDTCH
GNPSHYDKA VRRCCYRCPM GLFPTQCPQ
RPTDCRKQCE PDYYLDEADR CTACVTCSR
DLVEKTPCAW NSSRVCECRP GMFCSTSAVN
SCARCFHFSV CPAGMIVKFP GTAQKNTVCE
PASPGVSPAC ASPENCKEPS SGTIPQAKPT
PVSPATSSAS TMPVRGGTRL AQEAASKLTR
APDSPSSVGR PSSDPGLSPT QPCPEGSGDC
RKQCEPDYYL DEAGRCTACV SCSRDDLVEK
TPCAWNSSRT CECRPGMICA TSATNSCARC
VPYPICAAET VTKPQDMAEK DTTFEAPPLG
TQPDCNPTPE NGEAPASTSP TQSLLVDSQA
SKTLPIPTSA PVALSSTGKP
    
```

预测分子量	80 kDa
氨基酸	1 to 380
额外的序列信息	Fused to the Fc portion of human IgG1

技术指标

Our **Abpromise guarantee** covers the use of **ab157290** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Functional Studies
----	--------------------

	SDS-PAGE
形式	Lyophilized
补充说明	Binds Human CD30L (CD153).
制备和贮存	
稳定性和存储	Shipped at 4°C. Store at -20°C. Constituent: 99% PBS This product is an active protein and may elicit a biological response in vivo, handle with caution.
复溶	Reconstitute with 50µl sterile water to a concentration of 1 mg/ml. Further dilutions should be made with medium containing 5% fetal calf serum or a carrier protein. After reconstitution, prepare aliquots and store at -20°C. Avoid freeze/thaw cycles.
常规信息	
功能	Receptor for TNFSF8/CD30L. May play a role in the regulation of cellular growth and transformation of activated lymphoblasts. Regulates gene expression through activation of NF-kappa-B.
序列相似性	Contains 6 TNFR-Cys repeats.
翻译后修饰	Phosphorylated on serine and tyrosine residues.
细胞定位	Cytoplasm and Cell membrane.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors