

Recombinant human CD70 protein (His tag) ab222430

2 图像

描述	
产品名称	重组人CD70蛋白(His tag)
生物活性	Measured by its binding ability in a functional ELISA with Human CD27/TNFRSF7. The ED <sub>50</sub> range ≤ 0.5 ug/ml.
纯度	> 90 % SDS-PAGE. Affinity purified
内毒素水平	< 1.000 Eu/μg
表达系统	Baculovirus infected insect cells
Accession	<b>P32970</b>
蛋白长度	Protein fragment
无动物成分	No
性质	Recombinant
种属	Human
序列	ADPQRFAQAQQQLPLESLGWDVAELQLNHTGPQQDPRLYWQG GPALGRSF LHGPELDKGQLRIHRDGIYMVHIQVTLAICSSTTASRHHPTT LAVGICSP ASRSISLLRLSFHQGCTIASQRLTPLARGDTLCTNLGTLLP SRNTDETF FGVQWVRPHHHHHH
预测分子量	18 kDa including tags
氨基酸	39 to 193
标签	His tag C-Terminus
额外的序列信息	Extracellular domain (NP_001243).

技术指标	
Our <b>Abpromise guarantee</b> covers the use of <b>ab222430</b> in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	SDS-PAGE
形式	Liquid


制备和贮存

稳定性和存储	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.  pH: 7.40 Constituents: 10% Glycerol (glycerin, glycerine), PBS  This product is an active protein and may elicit a biological response in vivo, handle with caution.
--------	---


常规信息

功能	Cytokine that binds to CD27. Plays a role in T-cell activation. Induces the proliferation of costimulated T-cells and enhances the generation of cytolytic T-cells.
序列相似性	Belongs to the tumor necrosis factor family.
细胞定位	Membrane.

图片



Human CD27 Ligand/TNFRSF7 is coated at 10 ug/ml (100 µl/well) can bind Human CD27/TNFRSF7. The ED50 range ≤ 0.5 µg/ml.



15% SDS-PAGE analysis of 3 µg ab222430.

Lane 1: Non-reducing conditions.

Lane 2: Reducing conditions.

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors